

Surgical Gowns

Trusted cost-efficient solutions for your operating theatres



Surgical Gowns

Surgical gowns play a vital role in safeguarding both clinicians and patients within the operating theatre environment. Opting for single-use surgical gowns helps to minimise the risk infection.



Protection & Comfort

365 Healthcare surgical gowns are constructed from 5 layers of material to form SSMMS. The two top layers are spunbond with 2 middle layers of meltblown material and a final layer of spunbond material.

The material created is therefore strong, providing an efficient fluid and particulate barrier. For the user it is then efficient whilst enabling the fabric to breathe, providing comfort and freedom of movement.



Choice

The 365 Healthcare surgical gown range offers distinct areas of protection and reinforcement.

We understand the importance of protecting operating theatre teams and aim to match their needs. Our extensive selection of surgical gowns is designed to accommodate clinical preference in protection, safety, size and comfort. Our range caters for all team members with multiple sizes enabling comfort and freedom of movement during operations.


Surgical gowns from the 365 Healthcare range are offered in both lite and standard protection levels, plus reinforced and impervious ranges where more protection is needed. By offering a full variety of options your surgical team should be well catered for in the operating theatre.

Material Construction

Two layers of spunbond material

Two layers of meltblown material

One layer of spunbond material



Reassurance

The 365 Healthcare surgical gown range is fully compliant with EN13795.

By meeting EN13795 standards 365 Healthcare surgical gowns provide an effective barrier against contamination and infection for patients and clinicians.

Patient safety is at the centre of operating theatre objectives and by meeting the EN13795 standards the risk of surgical site infections is greatly reduced and leads to better outcomes for patients post surgery.


Quality is assured with the range as compliance to EN13795 means that healthcare professionals can use the surgical gown range knowing that it complies to a standard with specific performance requirements.

These rigorous standards encompass many test methodologies to ensure the surgical gowns perform in high-risk environments.

Benefits

- 1 Neck Hook & Loop Fastening**
For ease of adjustment and fit.
- 2 Ultrasonically Bonded Sleeves**
Maximum protection in critical areas.
- 3 Flexible Cuffs**
Reduces chance of skin debris being transferred to patient while maintaining dexterity.
- 4 Side Fastening**
Supports clinical protocols of gown donning.
- 5 Passcard**
Enables maintained sterility of the gown within the sterile area.
- 6 Anti-alcohol and anti-static treatment**
Enhances safety and clinical comfort.

All 365 Healthcare surgical gowns have colour coded printing on both the outer case and inner pouch insert, to easily identify the protection level of the surgical gown.



Specialised Protection

We provide surgical gowns for all operating theatre environments.

All our surgical gowns are constructed from 5 layers of SSMMS material for maximum strength and protection and meet the EN13795-1:2019 standard.

Each surgical gown is provided in a sterile pouch, wrapped in a sterile field complete with 2 generous hand towels. The gown is folded to support clinical protocols.


Each Pack Contains

1 x Sterile gown
1 x SMS wrap
2 x Hand towels

Standard Protection Gown

Used for moderate fluid procedures, ENT, Laparoscopic, Maternity (Labour Ward) and General Surgery.

Used in the majority of cases as a barrier to fluid. Suitable for moderate procedures with low blood or fluid spillage covering the majority of cases.




Size	NPC	Ref	Box Qty
Med	BWK85498	30-02S365	30
Lrg	BWK85486	30-03S365	30
XL	BWK85488	30-04S365	30
XL XLong	BWK85485	30-05S365	30
XXL	BWK85482	30-07S365	30
Med TL*	xxx	30-02S365/NC	30

* Thumb Loop

Standard Protection Lite Gown

Used for minimal low fluid procedures, Maxillofacial, Maternity (Labour Wards), Ophthalmic procedures, Cathlabs, Day Procedure Units and General Surgery.

Suitable for minimal fluid procedures whilst still providing a good barrier.



Size	NPC	Ref	Box Qty
Med	BWK85496	35-02LS365	30
Lrg	BWK85494	35-03LS365	30
XL	BWK85490	35-04LS365	30
XL XLong	BWK85492	35-05LS365	30


Zone Reinforced Protection Gown

Used for longer, higher fluid procedures. Major General Surgery, Orthopaedic & Urology.

Used where a procedure has a high level of fluid and thus offer the greatest level of protection. Appropriate for very wet, lengthy procedures where more fluid is expected such as Urology.

Size	NPC	Ref	Box Qty
Med	BWK85483	40-02ZR365	30
Lrg	BWK85412	40-03ZR365	30
Lrg*	BWK85499	40-23ZR365	30
XL XLong	BWK85500	40-04ZR365	30
XXL	BWK85497	40-05ZR365	30
XL*	BWK85495	40-24ZR365	30
XXL	BWK85493	40-07ZR365	30

* Extra long arthroscopy




Zone Impervious Protection Gown

Used for wet, lengthy procedures. Cardiovascular, Gynaecology, Major General, Orthopaedic.

Used where a procedure has a high level of fluid and thus offer the greatest level of protection. Appropriate for very wet, lengthy procedures where more fluid is expected such as Gynaecology & Orthopaedics.

Size	NPC	Ref	Box Qty
Med	xxx	50-02ZI365	30
Lrg	BWK85486	50-03ZI365	30
XL	xxx	50-04ZI365	30
XL XLong	xxx	50-05ZI365	30




Lite Zone Impervious Protection Gown

Used for moderate, wet, longer procedures. General Surgery, Neurosurgical, Orthopaedic, Cardiovascular and Maternity (C-Section).

Suitable for use where impervious protection is needed in critical areas. Cost-effective protection in operating theatre. Appropriate when moderate to high levels of fluid are expected, but the clinician requires a cooler gown.

Size	NPC	Ref	Box Qty
Med	BWK85491	90-02LZI365	30
Lrg	BWK85481	90-03LZI365	30
XL	BWK85489	90-04LZI365	30
XL XLong	BWK85487	90-05LZI365	30





How to Order

You can conveniently place orders for products via NHS Supply Chain or directly with us by calling **01530 278 111** or by emailing **365-sales@mediq.co.uk**

About Mediq

As part of Mediq UK, our Product Solutions division offers a diverse range of premium products, encompassing both our own brand range and partner offerings, tailored to meet your evolving needs. You can conveniently place orders for products directly via NHS Supply Chain or directly with us

providing flexible ordering which, alongside our excellent customer support, means you can gain greater efficiencies and avoid unnecessary and costly over supply and waste. We are compliant and used by healthcare institutions.

About the 365 Healthcare Brand

The 365 Healthcare brand, now part of the Mediq portfolio, is renowned for reducing the cost of care without clinical compromise.

At the heart of the brand's success is a customer-focused approach.

Each product within the 365 Healthcare range is developed through close collaboration with clinicians and healthcare providers. By engaging with hospitals, clinics, and private settings, the 365 Healthcare brand of products are designed to solve practical problems and fulfil real needs.

The reputation of the 365 Healthcare brand of products has been built upon both quality and affordability. The product ranges offer both tailored solutions and reliable supply from

our trusted manufacturing facilities, ensuring consistency and flexibility for our customers. With years of clinical use, our products are tried and tested in the market, providing the quality and support our customers expect.

As the 365 Healthcare brand of product expands, we are bringing compliant products into Europe while continuing to support our loyal UK customers. Our integration with Mediq, a leading European provider of medical devices and care solutions, strengthens our commitment to quality and innovation.

As a product brand within Mediq, 365 Healthcare remains dedicated to delivering cost-effective solutions without clinical compromise. Our products uphold the high standards established by both the 365 Healthcare brand and Mediq.

Additional Support Provided via:



Dedicated Account Managers and Clinical Specialists:
to support your unique requirements, preferences and challenges.



Implementation Excellence:
providing you different options from simple supply to full-service implementation. Full-service implementation encompasses audit of current products, product selection.



Savings Analysis:
Assisting you in savings reviews to demonstrate the most cost-effective solutions.



Dedicated Service Level Agreements:
The ability to provide bespoke SLAs.

Quick Reference Guide

Surgical Gown Selection

					
	Standard SSMMS	Standard Lite	Zone Reinforced	Zone Impervious	Lite Zone Impervious
Cardiovascular			●	●	●
ENT	●				
General Surgery	●	●	●	●	●
Gynaecology			●	●	
Laparoscopic	●				
Major General			●	●	
Maternity C-section					●
Maternity Labour Ward	●	●			
Maxillofacial	●	●			
Neurosurgical			●		●
Ophthalmic	●	●			
Orthopaedic			●	●	●
Patient Examination	●	●			
Urology			●	●	



at the heart of healthcare

Customer Services: 01530 278 111 | 365-sales@mediq.co.uk
www.mediq.co.uk

©2024 Mediq UK Limited. All rights reserved, 365 Healthcare® is a registered trademark of Mediq Ltd.
All other trademarks or brand names are the property of their respective owners.
Mediq UK Healthcare Ltd, Mediq House, Arundel Ave, Castle Donington, Derby DE74 2HJ.
Company Registration Number: 00062537, Reference Number: MED36501/0424