

365 Surgical Gowns



Surgical gowns are an essential part of protection in an Operating Theatre environment for both clinicians and patients. Using single use surgical gowns dramatically reduces the risk of infection.

365 Surgical Gowns are made from a polypropylene blend construction with 5 layers of melt and spun-blown material fibres. 365 Healthcare surgical gowns are fully compliant with EN13795. By meeting the EN13795 standards clinicians can be confident that they offer adequate protection in an Operating Theatre.

365 Surgical gowns provide:

- An efficient fluid and particulate barrier
- High level strength and absorbency
- Sound surgical bacterial barrier
- Lightweight, breathable fabric for comfort
- Ultrasonically melded sleeve seams
- Anti-Static and alcohol treated fabric
- Quality, cost effective gown solutions

This flyer contains information on 365 Surgical Gowns.

To view our full range of clothing and theatres products visit 365healthcare.com or scan the code below



Tel: 01905 778 365

Email: info@365healthcare.com



365 Gown Features

- Fully compliant with EN13795
- Hook and loop fastening at the neck to ensure secure and easily adjustable fit
- Soft, elasticated fabric cuffs for extra comfort
- Side tie-fastenings for accessibility
- Overlap sterile back, for hygiene and fit
- Colour coded print on outer case, to identify protection level
- Supplied with 2 generous hand towels
- Wrapped in a sterile field



365 Surgical Gown Range



REF	MATERIAL	SIZE	BOX QTY
-----	----------	------	---------

Standard Gowns - for low fluid procedures

35-02LS365	SSMMS Standard Lite	Medium	30
35-03LS365		Large	30
35-04LS365		Extra Large	30
35-05LS365		Extra Large, Extra Long	30
30-02S365	SSMMS	Medium	30
30-03S365		Large	30
30-04S365		Extra Large	30
30-05S365		Extra Large, Extra Long	30
30-07S365		Extra Extra Large	30
30-02S365/NC		Medium (thumb loop)	30

Reinforced Gowns - for moderate fluid procedures

40-02ZR365	SSMMS	Medium	30
40-03ZR365		Large	30
40-23ZR365		Large Extra Long Arthroscopy	30
40-04ZR365		Extra Large	30
40-05ZR365		Extra Large Extra Long	30
40-24ZR365		Extra Large Extra Long Arthroscopy	30
40-07ZR365		Extra Extra Large	30

Impervious Gown - for high fluid procedures

90-02LZI365	SSMMS Impervious Lite	Medium	30
90-03LZI365		Large	30
90-04LZI365		Extra Large	30
90-05LZI365		Extra Large Extra Long	30
50-02ZI365	SSMMS	Medium	30
50-03ZI365		Large	30
50-04ZI365		Extra Large	30
50-05ZI365		Extra Large Extra Long	30

Department	Standard Lite	Standard SMMS	Reinforced	Impervious
Cardiovascular				•
ENT		•		
General Surgery	•	•	•	•
Gynaecology				•
Laparoscopic		•		
Major General			•	•
Maxillofacial	•	•		
Neurosurgical				•
Ophthalmic	•	•		
Orthopaedic		•	•	•
Patient Examination	•	•		
Urology			•	•

Glossary of ISO 15223-1 Symbols

The following symbols are required by ISO 15223-1 to be included on packaging of medical devices. Relevant symbols are included above product tables throughout this brochure to aid you in your choice of an appropriate product.

Single Use
Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

CE Mark
CE marking is the medical device manufacturer's claim that a product meets the essential requirements of all relevant European Medical Device Directives.

REF
The manufacturer's catalogue number so that the medical device can be identified.

Manufacturer
The legal Manufacturer of the Medical Device

Date of Manufacture

Use-By Date

Batch Code
The manufacturer's batch code so that the batch or lot can be identified

Consult Instructions
Indicates there is additional information to what is included on the packaging which should be consulted prior to use

Sterile
Indicates a medical device that has been subjected to a sterilization process. The letter in the right box indicates the method of sterilisation, e.g. EO for ethylene oxide

Do Not Use if Package is Damaged Usually means if the packaging is wet or damaged it will no longer be sterile and/or effective