


Product Specifications Sheet

Specification Sheet No	SSN-PPETF-GWNS-010			
Product Name	Aprons / PPE Gowns			
Product Description	Patient Ward Gown with Hook & Loop Fastening			
Range	365 Healthcare			
Manufacturer	Mediq Healthcare UK Ltd			
Material	SMS, Nylon Hook and Loop			
Sterilisation	Non- Sterile			
Issue Date	13th September 2024	Version No	R01	

Intended use	A non-sterile, single use gown used for patients to wear during their stay in hospital.
--------------	---

Packaging							
Product code	Size	Colour	Length (cm) A	Width (cm) B	Length of Hook and Loop (cm) C	Items Per Pack	Qty – Shipper Carton
20-832365/NS	Medium	Dark blue	99±2	66±2	90±2	10	50
20-833365/NS	Large		105±2	71±2	90±2		
20-834365/NS	X Large		112±2	71±2	90±2		
20-835365/NS	XX Large		135±2	78±2	90±2		

Quality Certification	ISO 13485:2016		
Device Classification	PPE Category 1	Rule	N/A
Conformity Route	PPE (EU) 2016/425 with Minimal Risk		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	365 Healthcare Logo, Product name, Ref, Quantity, CE Mark, UKCA Mark, FSC, Lot, Non-sterile, Single use, Latex free, Store in cool dry place, Barcode, Manufacturing date, Expiry date, Manufacturer name & Address, ECREP name & Address, Importer name & Address		
Additional information	Products are Latex and phthalate free Products contain no hazardous substances Products contain no medicinal substances		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision, but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	13/09/2024	Update to Mediq Healthcare UK template from DSN698	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024