

GAMMEX[®] Latex

FORMULATION SHEET

PRODUCT DESCRIPTION: Latex surgical gloves - powder-free

FUNCTION	CHEMICAL
Base Material	Natural Rubber Latex
Antioxidant	Phenolic Antioxidant Benzimidazole Derivative
Accelerator	Xanthogen Polysulphide Sodium Dibutyldithiocarbamate Zinc Diethyldithiocarbamate
Vulcanizing agent	Sulphur
Activator	Zinc Oxide
Stabilizer	Potassium Hydroxide
Pigment	Titanium Dioxide

Confidentiality Notice: This document is confidential and contains proprietary information and intellectual property of Ansell Limited and/or its affiliates ("Ansell"). Neither this document nor any of the information contained herein may be reproduced or disclosed under any circumstances without the express written permission of Ansell. Please be aware that disclosure, copying, distribution or use of this document and the information contained therein is strictly prohibited.