



In application of the Regulation (EU) 2016/425 of 9 March 2016 concerning the harmonization of the Member States legislation relative to personal protective equipment, Centexbel Notified body 0493 authorized by the FPS Economy (Federal Public Services) has issued the following:

## EU TYPE EXAMINATION CERTIFICATE

**Nr. 032/2018/1898**

**This EU Type examination certificate is valid until 30 Oct 2023**

to: **Ansell Healthcare Europe nv, Brussel**  
for: **The Gloves Micro-Touch® Nitrile Accelerator-Free**

The personal protective equipment above mentioned satisfies the applicable essential safety requirements of the Regulation (EU) 2016/425.

For the argumentation, the following standards are used:

|                     |                                                                                                                                               |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| EN 420:2003+A1:2009 | Protective gloves – General requirements and test methods                                                                                     |
| EN 421:2010         | Protective gloves against ionizing radiation and radioactive contamination                                                                    |
| EN ISO 374-1:2016   | Protective gloves against dangerous chemicals and micro-organisms - Part 1: Terminology and performance requirements for chemical risks       |
| EN ISO 374-5:2016   | Protective gloves against dangerous chemicals and micro-organisms - Part 5: Terminology and performance requirements for micro-organisms risk |

If there is a former EC Type examination certificate according to the Directive 89/686/EEC this certificate remains valid until 21 April 2023 unless it expires before that date, for products that were manufactured before the issuance of this new EU Type examination certificate according to the Regulation (EU) 2016/425.

This is PPE of category III, subject to regular checks in accordance with article 19 of the European PPE Regulation. In agreement with the manufacturer's choice supervised product checks at random intervals shall be carried out to assess the quality of the final product (Module C2). The manufacturer must be able, on request, to present the test report of this quality control check. A first quality control check shall be performed at the latest on 31 Dec 2019 and at least be repeated with intervals of one year.

This declaration applies to the equipment as submitted in the type testing and described in the manufacturer's technical documentation (As described in 2016/425 Annex III) that is registered with number 10926.

Issued by Centexbel, Notified Body 0493, in Ghent, on 30 Oct 2018

Inge De Witte  
Certification Manager

Attached: 1 Annex




**CENTEXBEL • TEXTILE COMPETENCE CENTRE**

Technologiepark 7 • BE 9052 Gent • Belgium • phone +32 9 220 41 51 • fax +32 9 220 49 55 • [gent@centexbel.be](mailto:gent@centexbel.be) • [www.centexbel.be](http://www.centexbel.be)

VAT • BE 0459.218.289 • IBAN • BE 44 2100 4729 6545 • BIC • GEBABEBB

056-PROD

NOTIFICATION: EMBREUSE-BE 10050458-IND-1 VAN DE BELGISCHE WET VAN 21 JUNI 1987 (M.B. 1987, 1988, 1989, 1990, 1991, 1992, 1993, 1994, 1995, 1996, 1997, 1998, 1999, 2000, 2001, 2002, 2003, 2004, 2005, 2006, 2007, 2008, 2009, 2010, 2011, 2012, 2013, 2014, 2015, 2016, 2017, 2018, 2019, 2020, 2021, 2022, 2023) - EMBREUSE-BE 10050458-IND-1 VAN DE BELGISCHE WET VAN 21 JUNI 1987 (M.B. 1987, 1988, 1989, 1990, 1991, 1992, 1993, 1994, 1995, 1996, 1997, 1998, 1999, 2000, 2001, 2002, 2003, 2004, 2005, 2006, 2007, 2008, 2009, 2010, 2011, 2012, 2013, 2014, 2015, 2016, 2017, 2018, 2019, 2020, 2021, 2022, 2023)



## ANNEX

### EU TYPE EXAMINATION CERTIFICATE Nr. 032/2018/1898

#### 1. Applicant

Ansell Healthcare Europe nv  
Riverside Business Park Spey House  
Internationalelaan 55  
1070 Brussel  
Belgium

#### 2. Description

EN 421:2010



EN ISO 374-1 / Type C



EN ISO 374-5



VIRUS

#### 3. Materials and accessories

##### Gloves

- Micro-Touch® Nitrile Accelerator-Free



#### 4. Technical documentation

##### Summary test results

EN 420:2003+A1:2009 Gloves **Micro-Touch® Nitrile Accelerator-Free**

| Method           | Description                 | Result                  | Class   |
|------------------|-----------------------------|-------------------------|---------|
| EN 1413          | pH - textile                | PASS                    |         |
| EN 14362-1       | AZO dyes for colored gloves | PASS                    |         |
| EN 420 length    | Length                      | Fit for special purpose |         |
| EN 420 dexterity | Dexterity                   | PASS                    | Level 5 |

EN 421:2010 Gloves **Micro-Touch® Nitrile Accelerator-Free**

| Method                | Description                                           | Result                                    | Class   |
|-----------------------|-------------------------------------------------------|-------------------------------------------|---------|
| EN 420 dexterity      | Dexterity                                             | PASS                                      | Level 5 |
| EN 374-3              | Permeation K                                          | PASS                                      | Level 6 |
| EN 374-3              | Permeation T                                          | PASS                                      | Level 6 |
| EN 374-2              | Penetration                                           | PASS                                      |         |
| EN 388                | Mechanical requirements                               | Dexterity is the most important parameter |         |
| EN 421 4.3 method 5.1 | Attenuation efficiency and uniformity of distribution | /                                         | /       |

EN ISO 374-1:2016 Gloves **Micro-Touch® Nitrile Accelerator-Free**

| Method        | Description   | Result | Class   |
|---------------|---------------|--------|---------|
| EN 374-2      | Penetration   | PASS   |         |
| EN 16523-1    | Permeation K  | PASS   | Level 6 |
| EN 16523-1    | Permeation T  | PASS   | Level 6 |
| EN 374-4      | Degradation K | PASS   |         |
| EN 374-4      | Degradation T | PASS   |         |
| Type of glove |               | PASS   | Type C  |

EN ISO 374-5:2016 Gloves **Micro-Touch® Nitrile Accelerator-Free**

| Method            | Description                | Result | Class |
|-------------------|----------------------------|--------|-------|
| EN 374-2          | Penetration                | PASS   |       |
| ISO 16604 Proc. B | Protection against viruses | PASS   |       |

##### CENTEXBEL • TEXTILE COMPETENCE CENTRE

Technologiepark 7 • BE 9052 Gent • Belgium • phone +32 9 220 41 51 • fax +32 9 220 49 55 • gent@centexbel.be • www.centexbel.be

VAT • BE 0459.218.289 • IBAN • BE 44 2100 4729 6545 • BIC • GEBABEBB

## Description/Picture of article

Article      **The Gloves Micro-Touch® Nitrile Accelerator-Free**



The above picture is a general picture of the article. Possible variations of the above article can be present in the technical file.

**Note :**

Any modification in material, design, or other technical features must be brought to the attention of the Notified Body.