

Product Specification



Version: 01

Date of issue: 11/07/2022

Date of revision: 14/07/2022

Medipal Alcohol Skin Wipes with 2% Chlorhexidine

Medipal Alcohol Skin Wipes saturated with 70% IPA and 2% Chlorhexidine are intended for the removal of adhesive and plaster cast residue from the skin.



Product Description

Product Name	Medipal Alcohol Skin Wipe with 2% chlorhexidine (100 Single sachets)
Product Code	SSD682110MP
Languages	EN

Pack Type & Dimensions

Count	1 Wipe (per sachet)
Dimensions	100 x 100mm (wipe)
Weight	1.4g / sachet pack
Count per box	100 wipes (single sachets)
Colour of wipes	White
Wipes material description	18gsm Polypropylene

Solution

Solution Description	70% Isopropyl alcohol & 2% Chlorhexidine
Hazard identification	Flammable

Primary Packaging (Sachet)

Pouch Type	Aluminium / Paper laminate (50mm x 60mm)
Label	Printed pouch

Secondary Packaging (Box)

Carton Description	Printed cardboard box
Barcode	5025254 017678

Outer Packaging (Shipping Carton)

Case Description	Brown corrugated cardboard box (Double walled)
Barcode	5025254024591

Version: 01

Date of issue: 11/07/2022

Date of revision: 14/07/2022

Medipal Alcohol Skin Wipes with 2% Chlorhexidine

Packaging & Logistical Information

Product Name	Medipal Alcohol Skin Wipes with 2% Chlorhexidine
Product Code	SSD682110MP
Case Quantity	20 boxes per shipping carton
Case Dimensions [OD]	L430 x W290 x H160 mm
Case Weight	3.5kg
Case/Layer	8
Layers/Pallet	7
Cases/Pallet	56
Pallet Weight	216kg
Hazardous for Transport	Not applicable
Commodity Code for Export	3005909900

