

Product Specification

Version: 04:00

Date of issue: 22/02/2017

Date of revision: 03/08/2022

medipal

Medipal Alcohol Surface Wipes, with 2% chlorhexidine

Medipal Alcohol surface wipes with 2% chlorhexidine are intended for disinfecting hard surfaces of non-invasive medical devices within healthcare and care environments. For professional use only.

The wipes have been independently tested to a range of European testing standards (EN) which includes EN16615.



Product Description

Product Name	Alcohol Surface Wipe with 2% Chlorhexidine (200 Single Sachets)
Product Code	SSD684110MPCE
Languages	EN/FR/ES/AR

Pack Type & Dimensions

Count	1 wipe (per sachet)
Dimensions	195 x 105mm (wipe)
Weight	1.83g / sachet
Count per box	200 wipes (single sachets)
Colour of wipes	White
Wipes material description	18gsm polypropylene

Solution

Solution Description	70% Isopropyl alcohol & 2% Chlorhexidine gluconate
Hazard identification	Highly Flammable

Primary Packaging (Sachet)

Pouch Type	Aluminium / Paper (50mm x 55mm)
Label	Printed pouch

Secondary Packaging (Box)

Carton Description	Printed cardboard box
Barcode	5025254001899

Outer Packaging (Shipping Carton)

Case Description	Brown corrugated cardboard box (Double walled)
Barcode	5025254001912

Product Specification

Version: 04:00

Date of issue: 22/02/2017

Date of revision: 12/01/2021



Medipal Alcohol Surface Wipes, with 2% Chlorhexidine

Packaging & Logistical Information

Product Name	Alcohol Surface Wipe with 2% Chlorhexidine
Product Code	SSD684110MPCE
Case Quantity	10 boxes per shipping carton
Case Dimensions [OD]	L430 x W320 x H235 mm
Case Weight	5.3kg
Case/Layer	10
Layers/Pallet	4
Cases/Pallet	40
Pallet Weight	232kg
Hazardous for Transport	Not applicable
Commodity code for export	3005909900

