

Product Specifications Sheet

Specification Sheet No	SSN-PPETF-GWNS-009			
Product Name	Aprons / PPE Gowns			
Product Description	Premier Patient Examination Gown - Stockinette Cuff			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Polypropylene			
Sterilisation	Non- Sterile			
Issue Date	21 st Feb 2025	Version No	R03	

Intended use	To be used to protect the end user from patient care activities.
--------------	--

Packaging						
Product code	Colour	Gown Length (cm)	Gown Width (cm) circumference	Sleeve Length (cm)	Qty- Inner	Qty – Shipper Carton
5519	Blue	110 (±5)	138 (±1)	60 (±1.5)	1	50

Quality Certification	ISO 13485:2016		
Device Classification	PPE Category 1	Rule	N/A
Conformity Route	PPE (EU) 2016/425 with Minimal Risk		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	Premier Logo, Product name, Ref, Colour, Weight, Quantity, CE Mark, UKCA Mark, Non-sterile, Single use, Latex free, Store in cool dry place, Barcode, Lot, Manufacturing date, Expiry date, Manufacturer name & Address, ECREP name & Address		
Additional information	The product is long sleeve, stockinette cuff, neck and waist ties, one size, latex free, and non-sterile		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision, but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	06/06/2024	Update to Mediq Healthcare UK	Jack Wong
R02	07/06/2024	Update Mediq Healthcare UK address	Jack Wong
R03	21/02/2025	Update the product description	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.