

Product Specifications Sheet

Specification Sheet No	SSN-BAN(NS)-001			
Product Name	Bandages			
Product Description	Light Support Cotton Retention Bandages			
Range	Premier Brand			
Manufacturer	Mediq Healthcare UK Ltd			
Material	100% Cotton Yarn			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	Intended to be wrapped around a part of the body to secure a device such as a dressing and secured with tape or a safety pin.
--------------	---

Packaging						
Product code	Colour	Dimension (cm)	Ply	Items Per Pack	Qty- Inner	Qty – Shipper Carton
1966	Beige	5x4.5	12	12	144	1966
1967	Beige	7.5x4.5	12	12	144	1967
1968	Beige	10x4.5	12	12	144	1968
1969	Beige	15x4.5	6	12	72	1969

Quality Certification	ISO 13485		
Device Classification	I	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	None		
Markings	Premier Band logo, Ref, Quantity, Product name, Type 2 BP, Size, Barcode, Lot number, Date of Manufacture, Use by date, Manufacturer name & Address, CE Mark, Single use, Latex free, Medical device, Non-sterile, keep dry, keep away from sunlight, EC-REP		
Additional information	The product is latex free, in individual blister packaging.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	23/06/2023	Update to Mediq UK Brand replacing DSN 0177	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024