


Product Procedure Specifications Sheet

Specification Sheet No	SSN-FT1357-003		
Product Name	365-10404		
Product Description	Anaesthetic Pack		
Range	365 Healthcare		
Manufacturer	Shanghai ISO Medical Products Co Ltd (SIMP)		
Packaging	40pcs		
Sterilisation	EO Sterilisation		
Issue Date	02 nd May 2024	Version No	R01



Intended use	Surgical Dressing Kit is intended for the out-patient service clinical use and operation wound treatment.
--------------	---

Pack content				
Item	Product code	Quantity	Dimensions	Material
2 Part tray	01.03.01.0010	1	182X98X25MM	Polystyrene / Plastic
Drape Clear Fenestrated	03.08.00.0107	1	125X75CM with Fenestrated 12X7CM Adhesive Strip 16X11CM folded 15X19CM 35ply	Polyethylene / Plastic
Tissue Backed Sterile Field	01.08.00.0277	2	90x90cm	Polyethylene / Plastic
XRD Swab Gauze Green	03.01.00.0193	5	10x10cm 8ply	Cotton yarn with XRD thread
Syringe Luer Slip	01.00.00.0078	1	2..5ml	Piston Plug – Synthetic rubber Core bar and sheath – Polypropylene Lubricant – silicone oil Adhesive – Epoxy resin
Syringe Luer Lock	01.00.00.0294	1	3ml	
Syringe Luer Slip	01.00.00.0079	1	5ml	
Syringe Luer Slip	01.00.00.0081	1	10ml	
Syringe Luer Slip	01.00.00.0129	1	20ml	
Filter Needle	01.00.00.0517	1	18G x 4cm	Needle Stainless Steel Needle hub & shield - Polypropylene
Sponge Sticks Blue	01.00.00.0109	2		Polyurethane foam, PS
Hand Towels	03.08.00.0011	2	38x43cm	Pulp
Hypodermic Needle	01.00.00.0087	1	21g	Needle – Stainless Steel Needle hub & protective cap – Polypropylene
Hypodermic Needle	01.00.00.0152	1	23g	
Hypodermic Needle	01.00.00.0139	1	25g	

Quality Certification	ISO 13485		
Device Classification	Ila	Rule	N/A
Conformity Route	Medical Devices Directive (93/42/EEC)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Shelf life	5 Years if stored as recommended
Contraindications/Warnings	No Known contraindication
Markings	365 Healthcare, Ref, Product name, Lot Number, Quantity, Data Matrix, EC-REP, EO Sterile, Manufacturing date, Expiry date, CE Mark, Single use, Latex free, Medical device, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer Name & Address
Additional information	Latex free and sterile. Do not use if the packaging is broken.
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>	

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	02/05/2024	Update to Mediq UK Brand	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023