



Product Specifications Sheet

Specification Sheet No	SSN-FT1357-007			36590041  36590046 
Product Name	Wound Care Dressings			
Product Description	Reinforced Transparent IV Dressing			
Range	365 Healthcare			
Manufacturer	Shanghai ISO Medical Products Co., Ltd. (SIMP)			
Material	PU Film, Release paper, PET film, Solvent glue			
Sterilisation	EO Sterile			
Issue Date	3 rd June 2025	Version No	R03a	

Intended use	Securement Dressing
--------------	---------------------

Packaging								
Product code	Colour	A (mm)	B (mm)	C (mm)	D (mm)	Items per Pack	Qty – Inner	Qty – Shipper Carton
36590041	Transparent & white	93 ±3	65 ±3	60 ±2	51 ±2	1	100	2000
36590046		102 ±3	72 ±3	87 ±3	72 ±3	1	50	1000

Quality Certification	ISO 13485:2016		
Device Classification	Class IIa	Rule	7
Conformity Route	Medical Devices Directive (MDD 93/42/EEC)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare logo, Ref, Product name & dimension, Quantity, Barcodes, Lot, Manufacturing date, Expiry date, Single use, EO Sterile, CE Mark, Latex free, Medical device, consult instruction for use, keep dry, keep away from sunlight, Do not use if package is damaged, EC-REP, Manufacturer Name and Address, Distributor Name and Address		
Additional information	The product is sterile and latex free.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	29/05/2024	Update to Mediq Healthcare UK brand replacing DSN659	Elizabeth Ogunkolade
R02a	29/08/2024	Updated to remove product Transparent Ported IV Dressing 36590047	Elizabeth Ogunkolade
R03a	03/06/2025	Updated to revision 3	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.