


Product Specifications Sheet

Specification Sheet No	SSN-FT1357-010			 
Product Name	Wound Care			
Product Description	Transparent Island Dressings (film dressings with pads)			
Range	365 Healthcare			
Manufacturer	Shanghai ISO Medical Products Co Ltd (SIMP)			
Material	PU pad laminated with PE film, PE film thickness, PET tape Solvent adhesive & Release paper			
Sterilisation	EO Sterile			
Issue Date	29th August 2024	Version No	R01a	

Intended use	Film dressings with non-adherent pad indicated for low to moderate exuding acute, chronic, and surgical wounds, superficial cuts, burns and abrasions.
--------------	--

Packaging								
Product code	Colour	A(mm)	B(mm)	C(mm)	D(mm)	Items per Pack	Qty – Inner	Qty – Shipper Carton
36590019	Transparent & White	72 ±3	50 ±3	36 ±2	25 ±2	90	50	450
36590021		120 ±3	100 ±3	75 ±2	50 ±2	24		1200
36590022		155 ±3	85 ±3	116 ±3	38 ±2	12		600
36590023		200 ±3	100 ±3	150 ±3	53 ±2	12		600
36590024		250 ±3	100 ±3	200 ±6	55 ±2	12		600
36590025		300 ±3	100 ±3	250 ±6	55 ±2	12		600
36590026		95 ±3	85 ±3	68 ±2	38 ±2	24		1200
36590027		350 ±3	100 ±3	300 ±5	50 ±2	12		600

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare, Ref, Product name in 11 languages, Size, Image, Quantity, Barcodes, Lot, Manufacturing date, Expiry date, CE Mark, Single use, Latex free, EO Sterile, Medical device, Do not use if package is damaged, consult instruction for use, keep dry, keep away from sunlight, EC-REP, Manufacturer Name and Address, Distributor Name and Address		
Additional information	It is shower proof, and impenetrable to liquids, bacteria, and viruses to reduce the risk of wound infection. It is breathable allowing effective oxygen and moisture vapour exchange.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	29/08/2024	Update to Mediq Healthcare UK brand replacing DSN664	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024