


Product Specifications Sheet

Specification Sheet No	SSN-FT1357-011			
Product Name	Wound Care Dressing			
Product Description	Non-Adherent Dressing			
Range	365 Healthcare			
Manufacturer	Shanghai ISO Medical Products Co Ltd (SIMP)			
Material	100% Cotton and PE Film			
Sterilisation	EO Sterilisation			
Issue Date	29th August 2024	Version No	R01a	

Intended use	Low-adherent absorbent dressing used for the management of light to moderately excluding wounds including clean sutured wounds, abrasions, lacerations, and minor burns.
--------------	--

Packaging						
Product code	Colour	A mm	B mm	Items Per Pack	Qty- Inner	Qty – Shipper Carton
36519077	White	50±2.5	50±2.5	24	100	2400
36519078		100±5	100±5	18	50	900
36519079		100±5	200±10	18	25	450
36519080		200±10	300±10	6	25	150

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare, Ref, Product name, Size, Quantity, Barcodes, Lot, Manufacturing date, Expiry date, CE Mark, Single use, Latex free, Medical device, EO Sterile, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer Name and Address, EC-REP, Distributor Name and Address		
Additional information	Latex free and sterile. Do not use if the inner packaging is broken.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	29/08/2024	Update to Mediq Healthcare UK brand replacing DSN666	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.