

Product Procedure Specifications Sheet

Specification Sheet No	SSN-FT1357-014		
Product Name	Renal Dressing Pack		
Product Code	365-14116		
Range	365 Healthcare		
Manufacturer	SIMP (Shanghai ISO Medical Products Co., Ltd.)		
Packaging	110pcs per Outer		
Sterilisation	EO Sterile		
Issue Date	08 th April 2025	Version No	R02



Intended use	This range consists of general supplementary items which can be used on its own or in addition to the above procedure packs for the required clinical procedure.
--------------	--

Pack content				
Item	Product code	Quantity	Dimensions	Material
Clear Gallipot	01.03.00.0012	3	60ml	Polypropylene / Plastic
Non-woven Swab (4ply)	03.02.00.0001	2	5 x 5cm	70% rayon + 30% Polyester
Non-woven Swab (4ply)	03.02.00.0003	1	10 x 10cm	70% rayon + 30% Polyester
Cotton Sticks	01.00.00.0027	6	15cm	Stick – Wood, Tip - Cotton
2 part Tray (white)	01.03.01.0010	1	18.2 x 9.8 x 2.5cm	Polystyrene / Plastic
Green Forceps	01.00.00.0047	1	12.7cm	Polystyrene / Plastic
Dressing Towel	03.08.00.0011	1	40 x 38cm	Pulp 18gsm
Drapes	03.08.00.0028	2	50 x 50cm	20gsm PE + 25gsm Paper
Sterile Field	03.08.00.0027	1	50 x 50cm	20gsm PE + 25gsm Paper

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	Article 22
Conformity Route	Medical Device Directive (93/42/EEC As amended by 2007/47/EC) REGULATION (EU) 2023/607 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 Years – if stored as recommended above		
Contraindications/Warnings	No Known		
Markings	365 Healthcare logo, product reference, product name, Pack contents, Quantity, manufacture date, expiry date, lot number, barcode, single use only symbol, latex free symbol, medical device symbol, sterilised by ethylene oxide symbol, do not re-sterilise, do not use if packaging is damaged symbol, CE mark, keep dry symbol, keep away from sunlight symbol, manufacture address, EC rep details, Distributor address		
Additional information	This product is latex free. Do not use if the packaging is damaged. Product sterilised with ethylene oxide.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd.
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	03/10/2024	1 st version	Jack Wong
R02	08/04/2025	Update packaging QTY	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024