


Product Specifications Sheet

Specification Sheet No	SSN-FT1358-007		
Product Name	Surgical Gowns		
Product Description	Zone Reinforced Protection Surgical Gowns		
Range	365 Healthcare		
Manufacturer	Sino Protection (Hefei) Medical Products Co., Ltd		
Material	Gown: Blue AR SMMMMS Wrap: Blue SMS Reinforcement: White breathable PP + PE Towel: Polyester scrim		
Sterilisation	EO Sterile		
Issue Date	08 th April 2025	Version No	R02



A = Neck Circumference
B = Sleeve Length
E = Gown Circumference
G = Gown Length

Intended use	Intended to be used to prevent cross contamination between patients and end user
--------------	--

Packaging							
Product code	Size	A (cm) +2/-0	B (cm) +0/-2	E (cm) ±2	G (cm) ±2	Qty- pc/bag	Qty – pc/box
40-02ZR365	Medium	65	54	145	115	1 x Sterile Gown – Colour Blue, 1 x SMS Wrap, 2 x Hand towel	30
40-03ZR365	Large	65	58	156	122		
40-04ZR365	Extra Large	65	62	174	130		
40-05ZR365	Extra Large, Extra Long	65	65	176	151		
40-08ZR365	Extra, Extra Large, Extra Long	65	58	193	160		
40-09ZR365	Extra, Extra Large, Long	67	60	193	145		

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	Medical Devices Directive 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	No Known		
Markings	365 Logo, Ref, Size, Product name, Product list, Image, Quantity, EN 13795-1:2019, Zone reinforced protection: ZR, Barcode, NAFDAC number, Lot, Manufacture date, Expiry date, CE Mark, Single use, Latex free, Medical device, Sterile EO, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer name & Address, EC-REP		
Additional information	This product conforms to standard EN 13795-1:2019 and is latex free.		

Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	04/10/2024	1 st version	Jack Wong
R02	08/04/2025	Update dimension and add 40-09ZR365	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024