

Product Specifications Sheet

Specification Sheet No	SSN-FT1358-010		
Product Name	Surgical Gowns		
Product Description	Zone Impervious Protection Surgical Gown		
Range	365 Healthcare		
Manufacturer	Sino Protection (Hefei) Medical Products Co., Ltd		
Material	Gown: Blue AS + AR SMMMMS Wrap: Blue SMS Reinforcement: SPP+PE laminated Towel: Polyester scrim		
Sterilisation	EO Sterile		
Issue Date	08 th April 2025	Version No	R02

A = Neck Circumference
B = Sleeve Length
E = Gown Circumference
G = Gown Length

Intended use	Intended to be used to prevent cross contamination between patients and end user
--------------	--

Packaging							
Product code	Size	A (cm) +2/-0	B (cm) +0/-2	E (cm) ±2	G (cm) ±2	Qty- pc/bag	Qty – pc/box
50-02ZI365	Medium	65	54	145	115	1 x Sterile Gown – Colour Blue, 1 x SMS Wrap, 2 x Hand towel	30
50-03ZI365	Large		58	156	122		
50-04ZI365	Extra Large		62	174	130		
50-05ZI365	Extra Large, Extra long		65	176	151		

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	Medical Devices Directive 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	No Known		
Markings	365 Logo, Ref, Size, Product name, Product list, Image, Quantity, EN 13795-1:2019, Zone impervious protection: ZI, Barcode, Lot, Manufacture date, Expiry date, CE Mark, Single use, Latex free, Medical device, Sterile EO, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer name & Address, EC-REP		
Additional information	This product conforms to standard EN 13795-1:2019 and is latex free.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	04/10/2024	1 st version	Jack Wong
R02	08/04/2025	Update dimension	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024