

Product Specifications Sheet

Specification Sheet No	SSN-FT1390-001		
Product Name	Misc / Supplementary Items		
Product Description	30ml Measuring Cup		
Range	365 Healthcare		
Manufacturer	Finess Healthcare Group Co., Ltd		
Material	Polypropylene		
Sterilisation	Non-Sterile		
Issue Date	21st May 2024	Version No	R02



Intended use	A single use, disposable container suitable for use to dispense patient medication. Cup can be used for tablets and fluid which also has graduated marking.
--------------	---

Packaging						
Product code	Colour	Height $\pm$ 0.5 (mm)	Width $\pm$ 0.3 (mm)	Scale deviation (ml)	Items Per Pack	Qty – Shipper Carton
UN503	Transparent	41.5	39	<1.5	80pcs	4800pcs

Quality Certification	ISO 13485		
Device Classification	Class Im	Rule	1
Conformity Route	Medical Device Directive MDD 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	None		
Markings	365 Healthcare Logo, Ref Code, Quantity, Product name, Barcode, Lot, Manufacturing & Expiry date, CE Mark, Single use, Latex free, Non-Sterile, Medical device, keep dry, keep away from sunlight, do not use if box is damaged, Manufacturer name and address, Importer or Distributor's name and address, EC-REP		
Additional information	Product is latex free.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	15/12/2023	Update to Mediq UK Brand replacing DSN590 (product code UNH7-8005)	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.