


Product Specifications Sheet

Specification Sheet No	SSN-FT1390-002			
Product Name	Medicine Cup			
Product Description	Medicine Measure 60ml			
Range	365 Healthcare			
Manufacturer	Finess Healthcare Group Co., Ltd			
Material	Polypropylene			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	Intended to be used to for administering medicine to patient with graduated markings. Cup maybe used for tablets and fluid.
--------------	---

Packaging					
Product code	Colour	Weight (g) ±0.2	Size (mm) ±0.5	Qty- pc/bag	Qty – pc/box
UN504	Transparent	2.4	64.5 x 38	50	1000

Quality Certification	ISO 13485:2016		
Device Classification	Class I with measuring function	Rule	1
Conformity Route	Medical Devices Directives 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	3 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare Logo, Ref, Product name, Size, Quantity, Barcode, Lot, Date of Manufacture, Expiry date, CE Mark, Single use, Latex free, Non-Sterile, Medical device, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer name & Address, EC-REP, Distributor address		
Additional information	The product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	23/11/2023	New product specification sheet created	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.