


Product Procedure Specifications Sheet

Specification Sheet No	SSN-FT1390-015			
Product Name	Dressing Pack with Forceps			
Product Description	36510405			
Range	365 Healthcare			
Manufacturer	Wujiang Kangxin Surgical Dressings Co., Ltd.			
Packaging	120 Pieces			
Sterilisation	EO Sterile			
Issue Date	7th November 2024	Version No	R01a	

Intended use	Intended to be used to assist in minor Woundcare procedure
--------------	--

Pack content				
Item	Product code	Quantity	Dimensions	Material
Blue tissue backed drape	8802	1	50 x 50 cm	PE film and absorbent tissue
Blue tissue backed drape folded 24 ply	8702	1	50 x 50 cm	PE film and absorbent tissue
Yellow waste bag	2001	1	35 x 45 cm	High density polyethylene
White nonwoven balls (Medium)	6702	5	17 x 17 cm	70% viscose, 30% terylene
White nonwoven swab 4 ply	3102	5	7.5 x 7.5 cm	Viscose & terylene
Green plastic forceps	2303	2	12.5 cm	High impact polystyrene sheet
White 3-part tray	3604	1	18 x 9.7 x 2.3 cm	Polystyrene

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	Article 22
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored under recommended storage conditions		
Contraindications/Warnings	None		
Markings	Product Name, REF, 365 Logo, Lot Number, Manufacturing Date, Expiry Date, Medical Device, Single Use, EO Sterile, Single sterile barrier, do not use if package damaged, Keep dry and away from sunlight, CE, ECREP, Manufacturer name and address, Distributor name and address		
Additional information	Products are latex and phthalate free. Products contain no hazardous substances. Products contain no medical substances.		

Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	07/11/2024	Update to Mediq Healthcare UK brand replacing DSN752	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004-03	Procedure Pack Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024