

Product Specifications Sheet

Specification Sheet No	SSN-FT1580-011			 36590018  36590020
Product Name	Wound Care Dressings			
Product Description	Transparent Ported IV Dressing			
Range	365 Healthcare			
Manufacturer	Zhende Medical Co., Ltd			
Material	PU Film, Release paper, PET film, Solvent glue			
Sterilisation	EO Sterile			
Issue Date	24/07/2025	Version No	R02	

Intended use	Securement Dressing
--------------	---------------------

Packaging					
Product code	Colour	Product size (cm) ±3	Items per Pack	Qty – Inner	Qty – Shipper Carton
36590018	Transparent & white	7 x 9	50	40	2000
36590020		10 x 12	50	24	1200

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Devices Directive (MDD 93/42/EEC)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare logo, Ref, Product name in 11 other languages, Dimension, Quantity, Barcodes, Lot, Manufacturing date, Expiry date, Single use, EO Sterile, CE Mark, Latex free, Medical device, Do not re-sterilise, keep dry, keep away from sunlight, Do not use if package is damaged, EC-REP, Manufacturer Name and Address, Distributor Name and Address		
Additional information	The product is sterile and latex free. Moisture Vapour Transmission Rate >500g/m ² / 24hr		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	17/04/2024	Update to Mediq Healthcare UK brand replacing DSN659	Racheal Begg
R02	24/07/2025	Update with new picture	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.