


Product Specifications Sheet

Specification Sheet No	SSN-FT1580-013		
Product Name	Scrub Drapes		
Product Description	82-001365-Hip Set 1 82-002365-Hip Set 2 with Split Sheet 82-003365-Hip Set 3 with Split Sheet & Mayo		
Range	365 Healthcare		
Manufacturer	Xuchang Zhengde Environstar Medical Products Co., Ltd.		
Material	3-ply Coextruded material (Blue hydrophilic SPP + PE + white waterproof SPP), Antistatic SMS Blue, SS hydrophilic blue, Coextruded material (PE + white waterproof SPP), Blue OP-tape material, Dark blue SS reinforcement with PET material, Blue & grey PE cast film double corona thickness		
Sterilisation	EO Sterile		
Issue Date	28 th July 2025	Version No	R03



The diagram illustrates the dimensions for various medical drapes and covers. It includes a Top Drape (180 x 300 cm), a Bottom Drape (180 x 175 cm), a Split Drape (200 x 260 cm), a Utility Drape (150 x 175 cm), a Tape Strip (10 x 50 cm), a Trolley Cover (140 x 190 cm), a Mayo Cover (80 x 145 cm), a Hand Towel (39.37 x 48 cm), a Stockinette (30 x 120 cm), and a Utility Drape (75 x 90 cm).

Intended use	To be used to drape the patient to isolate a site from other areas of the patient's body or to contribute to reducing the risk of cross contamination.
--------------	--

Packaging							
Product code	Colour	Item	Dimension (cm)	Ply	Items per Pack	Qty – Inner	Qty – Shipper Carton
82-001365	Blue	1 x Top Drape	180 x 300 (±3)	3	1	4	8
		1x Bottom Drape	180 x 175 (±2)	3			
		1x Split Drape	200 x 260 (±3)	/			
		1x Tape Strip	10 x 50 (±1)	/			
		1x Trolley Cover	140 x 190 (±2)	/			
	White	4x Hand Towels	39.37 x 48 (±1)	2			
82-002365	Blue	1x Utility Drape	150 x 175 (±2)	3	1	3	6
		1x Top Drape	180 x 300 (±3)	3			
		1x Plastic Split Sheet	150 x 200 (±2)	/			
		1x Split Drape	200 x 260 (±3)	/			
		1x Tape Strip	10 x 50 (±1)	/			
		1x Trolley Cover	140 x 190 (±2)	/			
	White	4x Hand Towel	30.5 x 48 (±1)	2			
82-003365	Blue	1x Utility Drape	75 x 90 (±1)	3	1	4	8
		1x Top Drape	180 x 300 (±3)	3			
		1x Split Drape	200 x 260 (±3)	/			
		1x Mayo Cover	80 (+1/-0) x 145 (±2)	/			
		1x Stockinette Large	30 x 120 (±2)	/			
		2x Tape Strip	10 x 50 (±1)	/			
		2x Trolley Cover	140 x 190 (±2)	/			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Packaging							
Product code	Colour	Item	Dimension (cm)	Ply	Items per Pack	Qty – Inner	Qty – Shipper Carton
	White	4x Hand Towels	39.37 x 48 (±1)	2			

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	EU MDR 2017/745		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Healthcare, Ref, Product name in 11 languages, Dimension, Images, Quantity Barcodes, Lot, Manufacture date, Expiry date, CE Mark, Single use, Latex free, Medical device, EO Sterile, Do not use if package is damaged, keep dry, keep away from sunlight, EC-REP, Manufacturer Name & Address, Distributor Name & Address		
Additional information	This product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01a	08/05/2024	Update to Mediq Healthcare UK brand replacing DSN730	Elizabeth Ogunkolade
R02	24/07/2025	Update of dimension with tolerance	Jack Wong
R03	28/07/2025	Update the component of 82-001365	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023