


Product Specifications Sheet

Specification Sheet No	SSN-FT1580-022			
Product Name	Wound Care Dressings			
Product Description	Premierpore VP Non-woven Dressings			
Range	Premier			
Manufacturer	Zhende Medical Co.,Ltd.			
Material	PU coated non-woven, Acrylic glue, PE film, Viscose + polyester, Glassine paper			
Sterilisation	EO Sterile			
Issue Date	22 nd Jan 2025	Version No	R01	

Intended use	Film dressing with non-adherent pad indicated for low to moderate exuding acute, chronic, and surgical wounds, superficial cuts, burns and abrasions.
--------------	---

Packaging								
Product code	Colour	Dimension (cm)				Items Per Pack	Qty- Inner	Qty – Shipper Carton
		A	B	C	D			
5110N	White & Transparent	7	5	4	2.5	12	50	600
5111N		10	10	5	5	12	50	600
5112N		15	10	10	5	12	50	600
5113N		20	10	15	5	6	50	300
5114N		25	10	20	5	6	50	300
5115N		30	10	25	5	6	50	300

Quality Certification	ISO 13485:2016		
Device Classification	Class IS	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	Premier logo, Ref, Quantity, Product name, Size, Latex free, CE Mark, Single use, Medical device, EO Sterile, Do not use if package is damaged, Instruction for use, keep dry, keep away from sunlight, Barcodes, Lot, Manufacturing & Expiry date, Manufacturer Name and Address		
Additional information	Product is sterile, vapour permeable and Latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	22/01/2025	1 st version	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.