

Product Specifications Sheet

Specification Sheet No	SSN-FT1580-027		
Product Name	Surgical Drapes		
Product Description	Hip Drape with Side pouches-Elliptical		
Range	365 Healthcare		
Manufacturer	Xuchang Zhengde Environstar Medical Products Co.,Ltd.		
Material	2-ply Coextruded material dark blue hydrophilic, PE Film, Blue & Gray PE cast film double side corona, Transparent PEVA antistatic material, SMS Material		
Sterilisation	EO Sterile		
Issue Date	24 th July 2025	Version No	R01

The technical drawing illustrates the components and dimensions of the Hip Drape assembly. The main assembly, labeled 'Reinforced Hip Drape 300/230 x 354cm', consists of a large rectangular drape with a central opening for the hip. The overall dimensions are 300 cm in width and 354 cm in height. The central opening is 40 cm wide and 4.5 cm high. The side pouches are 230 cm wide. A 'Split Sheet 150x200 cm' is shown as a separate component, with dimensions 150 cm in width and 200 cm in height. A 'Hand Towel' is shown with dimensions 50 cm in width and 10 cm in height. A 'Trolley Cover' is shown with dimensions 190 cm in width and 140 cm in height.

Intended use	To be used to drape the patient to isolate a site from other areas of the patient's body or to contribute to reducing the risk of cross contamination.
--------------	--

Packaging							
Product code	Colour	Item	Dimension (cm)	Ply	Items per Pack	Qty – Inner	Qty – Shipper Carton
82-004365	Blue	Hip Drape	300/230 x 354 (±3)	2	1	4	8
		Plastic Split Sheet	150 x 200 (±2)	/			
		Tape Strip	10 x 50 (±1)	/			
		Trolley Cover	140 x 190 (±2)	/			
	White	4 x Hand Towels	39.37 x 48 (±1)	2			

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	No known		
Markings	365 Healthcare, Ref, Quantity, Barcodes, Lot, Manufacture date, Expiry date, CE Mark, Single use, Latex free, Medical device, EO Sterile, Do not use if package is damaged, keep dry, keep away from sunlight, EC-REP, Manufacturer Name & Address		
Additional information	This product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	24/07/2025	1sr version	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.