


## Product Specifications Sheet

|                        |                                                                      |            |     |                                                                                    |
|------------------------|----------------------------------------------------------------------|------------|-----|------------------------------------------------------------------------------------|
| Specification Sheet No | SSN-JAL007-001                                                       |            |     |  |
| Product Name           | Woundcare / Dressings                                                |            |     |                                                                                    |
| Product Description    | Non-Woven Island Dressing                                            |            |     |                                                                                    |
| Range                  | 365 Healthcare                                                       |            |     |                                                                                    |
| Manufacturer           | Allmed Medical Products                                              |            |     |                                                                                    |
| Material               | Non-Woven Fabric with 100% Polyester Release Paper, Solvent Adhesive |            |     |                                                                                    |
| Sterilisation          | EO Sterile                                                           |            |     |                                                                                    |
| Issue Date             | 17 <sup>th</sup> April 2024                                          | Version No | 01a |                                                                                    |

|              |                                                                      |
|--------------|----------------------------------------------------------------------|
| Intended use | Adhesive Island dressing for the treatment of low to moderate wounds |
|--------------|----------------------------------------------------------------------|

| Packaging    |        |                               |            |                |            |                      |
|--------------|--------|-------------------------------|------------|----------------|------------|----------------------|
| Product code | Colour | Dimension of product (±0.3cm) |            | Items Per Pack | Qty- Inner | Qty – Shipper Carton |
|              |        | Length (cm)                   | Width (cm) |                |            |                      |
| 36590001     | White  | 5                             | 7.2        | 1              | 50         | 3000                 |
| 36590002     |        | 10                            | 10         |                | 50         | 1200                 |
| 36590003     |        | 10                            | 15         |                | 50         | 600                  |
| 36590004     |        | 10                            | 20         |                | 25         | 300                  |
| 36590005     |        | 10                            | 25         |                | 25         | 300                  |
| 36590016     |        | 8                             | 15         |                | 50         | 600                  |
| 36590017     |        | 8                             | 10         |                | 50         | 600                  |
| 36590039     |        | 6                             | 8.3        |                | 60         | 3600                 |

|                                                                                                                   |                                                                                                                                                                                                                                                                                                                              |      |   |
|-------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|---|
| Quality Certification                                                                                             | ISO 13485                                                                                                                                                                                                                                                                                                                    |      |   |
| Device Classification                                                                                             | Class Is                                                                                                                                                                                                                                                                                                                     | Rule | 4 |
| Conformity Route                                                                                                  | Medical Devices Regulation (EU) 2017/745<br>Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)                                                                                                                                                                                                      |      |   |
| Storage conditions                                                                                                | Store in a cool dry place. Keep away from heat, direct sunlight.                                                                                                                                                                                                                                                             |      |   |
| Shelf life                                                                                                        | 5 years with the recommended storage conditions.                                                                                                                                                                                                                                                                             |      |   |
| Contraindications/Warnings                                                                                        | None                                                                                                                                                                                                                                                                                                                         |      |   |
| Markings                                                                                                          | 365 Healthcare, Ref, Product name, Dimension, Quantity, Barcodes, Lot, CE Mark, Manufacture date, Expiry date, Single use, Latex free, Medical device, EO Sterile, Do not re-sterilise, Do not use if package is damaged, keep dry, keep away from sunlight, EC-REP, Manufacturer Name & Address, Distributor Name & Address |      |   |
| Additional information                                                                                            | Products are Latex free, contain no hazardous substances, MVTR>500g/24 hours                                                                                                                                                                                                                                                 |      |   |
| Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices. |                                                                                                                                                                                                                                                                                                                              |      |   |

| Document Reference | Document Name               | Reviewed By | Approved By  | Revision Number | Issue Date                  |
|--------------------|-----------------------------|-------------|--------------|-----------------|-----------------------------|
| DMS 04-004         | Product Specification Sheet | Jack Wong   | Racheal Begg | R02             | 03 <sup>rd</sup> April 2024 |

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

|                                                            |                                  |
|------------------------------------------------------------|----------------------------------|
| 365 Healthcare                                             | Mediq House                      |
| T: +44 (0) 1332 413 700                                    | Arundel Avenue, Castle Donington |
| W: <a href="https://mediq.co.uk/">https://mediq.co.uk/</a> | Derby, DE74 2HJ                  |

365 Healthcare is a trading division of Mediq Healthcare UK Ltd  
Registered in England No. 62537

## Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

| Revision Number | Revision Date | Reason for Revision                       | Updated by           |
|-----------------|---------------|-------------------------------------------|----------------------|
| R01a            | 17/04/2024    | Update to Mediq UK Brand replacing DSN658 | Elizabeth Ogunkolade |
|                 |               |                                           |                      |
|                 |               |                                           |                      |

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.