


Product Specifications Sheet

Specification Sheet No	SSN-JAL007-007		
Product Name	Wound Care Dressings		
Product Description	Wound Dressings Pads		
Range	Premier		
Manufacturer	Allmed Medical Products Co., Ltd		
Material	Repellent Web 100% polypropylene, Inner (core) Pulp 100% wood pulp		
Sterilisation	ETO Sterile		
Issue Date	22 nd Oct 2024	Version No	R01



The image shows a rectangular, translucent Premier Wound Care Dressing Pad. Above the pad is the Premier logo with the tagline 'Healthcare Excellence'. The pad is placed on a white background with a grid. Red double-headed arrows indicate the dimensions: 'Width' for the vertical measurement and 'Length' for the horizontal measurement. The pad has a textured, non-adhesive outer layer and a smoother, adhesive inner layer.

Intended use	A single use, multi-layered wound care dressing pad used to absorb wound exudates and bodily fluids.
--------------	--

Packaging						
Product code	Colour	Dimension (cm)	Avg. Weight per pad (g)	Items Per Pack	Qty- Inner	Qty – Shipper Carton
2002N	White	10x10	3.0	24	50	1200
2004N	White	20x10	7.5	24	25	600
2006N	White	20x20	15.5	12	25	300
2008N	White	40x20	28.5	12	15	180

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Device Directive 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	Premier Logo, Ref, Code, Quantity, Product name, Dimension, Colour, CE Mark, Single use, Latex free, medical device, EO Sterile, Do not use if package is damaged, keep dry, keep away from sunlight, GS1 Barcodes, Lot, Manufacturing & Expiry date, EC-REP, Manufacturer name and address		
Additional information	Latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	22/10/2024	1 st version	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.