

Product Specifications Sheet

Specification Sheet No	SSN-JAL007-008		
Product Name	Wound Care Dressings		
Product Description	Premierpore Non-woven Dressings		
Range	Premier		
Manufacturer	Allmed Medical Products		
Material	Release Paper, PE, Viscose and Polyester, Adhesive, Non-woven Fabric		
Sterilisation	EO Sterile		
Issue Date	22 nd Jan 2025	Version No	R01a



Intended use	Film dressing with non-adherent pad indicated for low to moderate exuding acute, chronic, and surgical wounds, superficial cuts, burns and abrasions.
--------------	---

Packaging					
Product code	Finished product Dimension L x W (cm)	Absorbent Pad Dimension L x W (cm)	Items Per Pack	Qty- Inner	Qty – Shipper Carton
5101N	5x7	2.5x4	100	16	1600
5102N	10x10	5x5	100	12	1200
5103N	10x15	5x10	50	12	600
5104N	10x20	5x15	50	9	450
5105N	10x25	5x20	50	9	450
5106N	10x30	5x25	50	9	450
5107N	10x35	5x30	50	9	450

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	Premier logo, Ref, Product name & dimension, Quantity, Barcodes, Lot, Manufacturing date, Expiry date, CE Mark, Single use, Latex free, Medical device, EO Sterile, Do not use if package is damaged, keep dry, keep away from sunlight, Manufacturer Name and Address, EC-REP		
Additional information	Product is sterile Latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	22/01/2025	1 st version	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.