

Product Specifications Sheet

Specification Sheet No	SSN-JFRU06-015			
Product Name	Suction			
Product Description	Non-Conductive Bubble Tube			
Range	Universal			
Manufacturer	Romsons International			
Material	Polyvinyl Chloride			
Sterilisation	Non-sterile			
Issue Date	21 st May 2025	Version No	R01	

Intended use	Intended for connecting a suction source to suction cannisters, suction catheters, Yankauer tubes and other suction devices.
--------------	--

Packaging					
Product Code	Colour	Length (Approx.)	Diameter (Approx.)	Items Per Pack	Qty – Shipper Carton
UN730	Transparent	50m	5mm	1	6 packs
UN740		30m	5mm	1	6 packs
UN750		50m	6mm	1	6 packs
UN760		30m	7mm	1	6 packs

Quality Certification	ISO 13485:2016		
Device Classification	Class I	Rule	1
Conformity Route	EU MDR 2017/745		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years from manufacturing date if stored as recommended		
Contraindications/Warnings	Do not Use if packing is damaged. Tear Open pack Discard after Single Use. For Single Use Only.		
Markings	Universal Logo, Ref, Product name, Lot, Barcodes, Manufacture date, Expiry date, CE Mark, Latex free, Non-sterile, DEHP free, Medical device, keep dry, keep away from sunlight, Manufacturer, Distributor		
Additional information	Product is Latex-free, conductive. Easy to cut through the bubble to create the required length with a tapered end to fit any suction catheter or suction tip.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	21/05/2025	New SSN Mediq Healthcare UK brand Template	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023