

Product Specifications Sheet

Specification Sheet No	SSN- JSA001-001			 UN42450  UN42460  UN42470
Product Name	Uniprotect Respirator			
Product Description	UN42450 Uniprotect Unvalved FFP2 Respirator UN42460 Uniprotect Unvalved FFP3 Respirator UN42470 Uniprotect valved FFP3 Respirator			
Range	Universal			
Manufacturer	San Huei United Co. Ltd			
Material	Outer layer: Polypropylene Support layer: Polypropylene Filter Media: High Filtration Efficiency Melt-Blown (Polypropylene) Inner Layer: Polypropylene Nose Piece: PE Encapsulated Iron Wire. Headband: Latex-free Synthetic Rubber. Valve Housing: Polypropylene Valve Membrane: Synthetic Rubber			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	FFP2 NR D respirators can provide protection against solid and liquid aerosols and smoke. FFP3 NR D respirators can provide protection against solid and liquid aerosols and smoke. Features exhalation valve for extra comfort in hot environment and reduce hot air build up. (Only available on UN42470)
--------------	--

Packaging				
Product code	Colour	Dimension (cm) (Approximately)	Items Per Pack	Qty – Shipper Carton
UN42450	White	Width 85mm (upper end) and 200mm (lower end) Height 87mm	20pcs	240pcs
UN42460	White		20pcs	240pcs
UN42470	White		20pcs	240pcs

Quality Certification	ISO 13485		
Device Classification	PPE III	Rule	N/A
Conformity Route	Personal Protective Equipment Regulations (EU) 2016/425		
Storage conditions	Store in a cool dry place away from heat and direct sunlight. Storage temperature between -10°C to + 50°C		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	None		
Markings	Universal Logo, Ref Code, Quantity, Product name, Lot, Manufacturing & Expiry date, CE Mark, Single use, Latex free, Non-Sterile, Keep dry, keep away from		

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

	sunlight, IFU, Storage temperature between -10°C to + 50°C, EN 149:2001 + A1:2009, Manufacturer name and address, Importer or Distributor's name and address
Additional information	Product is latex free. Please read Instruction for Use before use. This respirator is limited to single shift use only (with a maximum of 8 hours) and shall not be re-used.
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>	

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	12/12/2023	SSN for supplier CE and 365 branding	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024