

Product Specifications Sheet

Specification Sheet No	SSN-MASK(NS)-002			
Product Name	Facemasks & Respirators			
Product Description	FluidProtect Anti-fog Surgical Facemask with Ties (Type IIR) FluidProtect Anti-fog Surgical Facemask with ear loops (Type IIR) FluidProtect Anti-fog Surgical Facemask with Visor and Ties (Type IIR) FluidProtect Anti-fog Surgical Facemask with Visor, Ear Loops (Type IIR)			UN49205
Range	Universal			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Outer or observer layer: Fluid resistant Spunbond Polypropylene fibre (hydrophobic) 2nd layer: Fluid resistant perforated PE film (PPF) Filter or 3rd Layer: 100%PP Meltblown filtration fibre (hydrophobic) Inner Layer: Fluid resistant Spunbond Polypropylene fibre (hydrophobic) Ear loops (UN49210 & UN49250): Spandex polyurethane coated with nylon Ties (UN49205 & UN49245) : PP Spunbond (hydrophobic)			UN49210
Sterilisation	Non-sterile			
Issue Date	20 th Nov 2024	Version No	R03	
				UN49245
				
				UN49250

Intended use	Intended to be placed over the mouth and nose of a healthcare worker to permit normal breathing whilst protecting the patient from the transfer of particles and contamination during medical procedures and surgery. Type IIR mask are fluid resistant.
--------------	--

Packaging							
Product code	Remark	Colour	Mask Body Dimension (mm) L x W	Ear Loop/Ties length (mm)	Items per Pack	Qty – Inner	Qty – Shipper Carton
UN49205	With ties	Orange	175 x 95 (±2)	900±5	6	50	300
UN49210	Ear loops			165 ±5	6	50	300
UN49245	Visor, ties			900 ±5	4	25	100
UN49250	Visor, Ear loops			165 ±5	4	25	100

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	3 years if stored as recommended		
Contraindications/Warnings	None		

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Markings	Universal Logo, Product name, Ref, Quantity, Lot number, Manufacturing date, Expiry date, Barcode, Manufacturer name & Address, Non-sterile, CE Mark, Single use, Latex free, medical device, keep dry, keep away from sunlight, EC-REP
Additional information	Microbial Cleanliness ≤ 30 cfu/g, splash resistance pressure ≤ 16.0 kPa
<i>Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.</i>	

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	23/06/2023	Update to Mediq UK Brand replacing DSN 829	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	20/11/2024	Update the product life to 3 years according to supplier specification	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024