

Product Specifications Sheet

Specification Sheet No	SSN-PPETF-BED-001		
Product Name	Bedding set		
Product Description	Premier Waterproof Mattress Cover		
Range	Premier		
Manufacturer	Mediq Healthcare UK Ltd.		
Material	Embossed cast LLDPE / LDPE		
Sterilisation	Non-Sterile		
Issue Date	16 th July 2024	Version No	R01



Intended use	To be used to protect the mattress from liquids and bodily fluids
---------------------	---

Packaging							
Product code	Colour	Length (approx.)	Width (approx.)	Height (approx.)	Thickness (approx.)	Qty – Inner	Qty – Shipper Carton
2501	Blue	210cm	90cm	20cm	0.030mm	10pcs	100pcs

Quality Certification	ISO 13485		
Device Classification	PPE Category I (Protection from Minimal Risk)	Rule	N/A
Conformity Route	(EU) 2016/425 Minimal Risk PPE Category 1		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	3 years if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	Premier Logo, Ref, Quantity, Product Description, Dimension, Barcode, Lot Number, Manufacture and Expiry Date, CE Mark, UKCA, Non-Sterile, Single use, Latex Free, keep dry, keep away from sunlight, Manufacturer Address, EC-REP		
Additional information	The product is latex free and non-sterile		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	16/07/2024	Update to Mediq Healthcare UK brand replacing DSN036	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.