


Product Specifications Sheet

Specification Sheet No	SSN-PPETF-GWNS-002			
Product Name	Aprons			
Product Description	Premier Disposable & Plastic Aprons			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	LDPE			
Sterilisation	Non- Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	For use in areas requiring protection against low level hazards, the effect of which are superficial and minimal
--------------	--

Packaging								
Product code	Colour	Apron Dimension L x W (cm)	Strip Dimensions L x W (cm)	Neck Dimension L x W (cm)	Thickness (micron)	Items Per Pack	Qty- Inner	Qty – Shipper Carton
2818N	Transparent	117x69 (± 5%)	40x6 (± 5%)	21x15 (± 5%)	18 (± 5%)	5	200	1000
2820N	White							
2822N	Blue							
2824N	Green							
2826N	Red							
2828N	Yellow							
2830N	Pink	135x69 (± 5%)						
2860N	White							
2862N	Blue							

Quality Certification	ISO 13485:2016		
Device Classification	PPE Category 1	Rule	N/A
Conformity Route	PPE (EU) 2016/425 with Minimal Risk		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	Premier Logo, Product name, Ref, Colour, Dimension, Thickness, Quantity, CE Mark, UKCA Mark, Nonn-sterile, Single use, Latex free, Store in cool dry place, Barcode, Lot, Manufacturing date, Expiry date, Manufacturer name & Address, ECREP name & Address		
Additional information	The product is latex free and non-sterile		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	26/09/2023	Update to Mediq Healthcare UK brand replacing DSN	Jack Wong
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024