


Product Specifications Sheet

Specification Sheet No	SSN-PPETF-GWNS-007		
Product Name	Aprons		
Product Description	Premier Heavy Duty Apron Green Roll		
Range	Premier		
Manufacturer	Mediq Healthcare UK Ltd		
Material	LDPE		
Sterilisation	Non- Sterile		
Issue Date	08 May 2024	Version No	R01



The diagram illustrates the dimensions of a green apron. On the left, a person is shown wearing the apron, with a vertical dimension line labeled 'A' indicating the 'Apron length'. To the right, a detailed line drawing of the apron shows various dimensions: 'B' for 'Apron width' at the top, 'C' for 'Strip length' on the left side, 'D' for 'Neck loop length' in the center, 'E' for 'Neck loop width' on the right side, and 'F' for 'Apron width' at the bottom. A small box labeled 'Apron length' is also present near the person's feet.

Intended use	For use in areas requiring protection against low level hazards, the effect of which are superficial and minimal
--------------	--

Packaging								
Product code	Colour	Apron Dimension length x width (cm)	Strip Dimension length x width (cm)	Neck Dimension length x width (cm)	Thickness (micron)	Items Per Pack	Qty- Inner	Qty – Shipper Carton
2884	Green	135x69 (± 5%)	40x6 (± 5%)	21x15 (± 5%)	50 (± 5%)	10	75	750

Quality Certification	ISO 13485:2016		
Device Classification	PPE Category 1	Rule	N/A
Conformity Route	PPE (EU) 2016/425 with Minimal Risk		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	N/A		
Markings	Premier Logo, Product name, Ref, Colour, Dimension, Thickness, Quantity, CE Mark, UKCA Mark, Non-sterile, Single use, Latex free, Store in cool dry place, Barcode, Lot, Manufacturing date, Expiry date, Manufacturer name & Address, ECREP name & Address		
Additional information	The product is latex free and non-sterile		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision, but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	08/05/2024	Update to Mediq Healthcare UK brand	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.