

Product Specifications Sheet

Specification Sheet No	SSN-PPETF-MOBCAPS-001		
Product Name	Headwear		
Product Description	Premier Mob Cap		
Range	Premier		
Manufacturer	Mediq Healthcare UK Ltd.		
Material	Non-Woven Polypropylene		
Sterilisation	Non-Sterile		
Issue Date	19 th September 2025	Version No	R03



Intended use	Ideal for hygiene-based environments including healthcare, beauty, food prep and catering, industrial and cleanrooms.
--------------	---

Packaging						
Product code	Colour	Name	Dispenser	Dimension (± 2cm)	Qty - Inner	Qty – Shipper Carton
2906	White	MOB CAP – Bouffant Style	Dispenser Box	56	100	1000
2907	Green					
2908	Blue					
2907-M	Green	Surgical Cap – Bouffant Style		50	150	600
2907-L	Green			53	150	600
2907-XL	Green			60	150	600
2907SW-L	Green	Surgical Cap – Bouffant Style with Sweatband		53	120	480
3906	White	Premier Mob Cap – Bouffant Style	Dispenser Bag	56 ± 2	100	1000
3908	Blue					

Quality Certification	ISO 13485		
Device Classification	PPE Category I	Rule	N/A
Conformity Route	PPE Regulation (EU) 2016/425		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	No Known		
Markings	Premier logo, Ref, Quantity, Product name, Colour, One size, Barcode, Lot, Manufacture date, Expiry date, CE Mark, UKCA, Non-sterile, Single use, Latex free, keep dry, keep away from sunlight, Manufacturer name & Address, EC-REP		
Additional information	The product is latex free, Fully Elasticated Edge, and non-sterile.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	28/03/2024	Update to Mediq UK Brand replacing DSN067	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	19/08/2025	Update with 2907-M, 2907-L, 2907-XL, 2907SW-L	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024