

Product Specifications Sheet

Specification Sheet No	SSN-PPETF-OSH-001			
Product Name	PE Overshoes			
Product Description	Premier Disposable Blue Overshoes			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd.			
Material	Polyethene			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	For use in areas requiring protection against mechanical action, the effect of which are superficial and light industrial.
--------------	--

Packaging								
Product code	Size	EU Shoe Size	Dimension (Approximately) L x W (cm)	Colour	Thickness (Approximately) (mm)	Qty - Bundle	Qty - Polybag	Qty - Shipper Carton
2401N	Small	33-37	36 x 15	Blue	0.04	25	100	2000
2402N	Large	38-42	41 x 15					

Quality Certification	ISO 13485		
Device Classification	PPE Category I	Rule	N/A
Conformity Route	PPE Regulation (EU) 2016/425		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended above		
Contraindications/Warnings	None		
Markings	Premier logo, Ref, Quantity, Product name, Colour, Size, Barcode, Lot, Manufacture date, Expiry date, CE Mark, UKCA, Medical device, Non-sterile, Single use, Latex free, keep dry, keep away from sunlight, Manufacturer name & Address, EC-REP		
Additional information	Premier overshoes have elasticated tops for a snug and comfortable fit. They are anti-static and anti-slip and suitable for use in clean conditions. The product is latex free and non-sterile.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	04/12/2023	Update to Mediq UK Brand replacing DSN035	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.