

Product Specifications Sheet

Specification Sheet No	SSN-PSW(NS)-001			
Product Name	Plain swabs			
Product Description	White Non-woven Swabs			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	50% Viscose and 50% Polyester			
Sterilisation	Non-sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	Intended to be used for non-invasive procedures such as skin cleansing.
--------------	---

Packaging						
Product code	Colour	Dimension (cm) Wrap x Weft	Items Per Pack	Ply	Qty- Inner	Qty – Shipper Carton
1910	White	5X5 (±5%)	4	200	40	8000
1915		7.5X7.5 (±5%)	4	200	20	4000
1918		7.5X7.5	6	100	50	5000
1920		10X10 (±5%)	4	200	20	4000
1922		10X10	6	100	50	5000

Quality Certification	ISO 13485		
Device Classification	I	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years if stored as recommended above		
Contraindications/Warnings	No known		
Markings	Premier logo, code ref, product name, ply, dimension, quantity colour, non-sterile, CE mark, single use, latex free, medical device, keep dry, keep away from sunlight, barcodes, manufacturing date, expiry date, lot number, manufacturer name and address, EC-REP		
Additional information	The product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	23/06/2023	Update to Mediq UK Brand replacing DSN 014	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024