

Product Specifications Sheet

Specification Sheet No	SSN-SA16-002			
Product Name	Misc / Supplementary Items			
Product Description	Red Polythene Laundry Bag			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Medium density polyethylene (MDPE)			
Sterilisation	Non- Sterile			
Issue Date	9th May 2024	Version No	R01	

Intended use	Intended to provide a hygienic way of isolating and transporting dirty or soiled linen.
--------------	---

Packaging							
Product code	Colour	Length (mm) approx.	Width (mm) approx.	Gusset (mm) approx.	Thickness- (micron)	Qty- Inner	Qty – Shipper Carton
2845	Red	762	711	457	10 (±15%)	50	200

Quality Certification	ISO 13485:2016		
Device Classification	N/A	Rule	N/A
Conformity Route	N/A		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	2 years if stored as recommended above		
Contraindications/Warnings	This product is not suitable for tunnel washers and tumble-drying machines. This sack has a water-soluble strip susceptible to damp conditions. Please ensure the tops of the sacks are secured using the water-soluble string tie to be found in the gusset		
Markings	Premier Logo, Product name, Ref, Quantity, Non-sterile, Single use, Latex free, keep dry, keep away from direct sunlight, Barcode, Lot, Manufacturing date, Expiry date, Manufacturer name & Address,		
Additional information	The product is latex free and non-sterile. Load the bag directly into the washing machine to eliminate the need to open bags to sort and/or load the contaminated linen as well as the potential exposure to the handlers. During the agitation cycle, the water-soluble seam and tie closure dissolve, releasing the contents of the bag into the washing solution. At the conclusion of the wash cycle, the used and disinfected portion of the polyethylene bag is removed and discarded or recycled. Minimum operating temperature for soluble strip is 23°C.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision, but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Ashley Miller	Racheal Begg	R01	12 th June 2023

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	09/05/2024	Update to Mediq Healthcare UK brand replacing DSN047	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.