


Product Specifications Sheet

Specification Sheet No	SSN-SCAP(NS)-006		
Product Name	Surgical Caps		
Product Description	Operating Caps – Tie Backs		
Range	Universal		
Manufacturer	Mediq Healthcare UK Ltd		
Material	Non-woven Polypropylene		
Sterilisation	Non Sterile		
Issue Date	21st May 2024	Version No	R02



Intended use	A non-sterile head covering designed as a cap intended to be worn by surgical staff during an operation to protect both the patient and themselves from the transfer of microorganisms, bodily fluids, and particulate material.
--------------	--

Packaging						
Product code	Colour	Dimensions (cm)			Qty- Inner	Qty – Shipper Carton
UN91250	Blue	A	Bottom Width	48 (±2)	100	600
		B	Tie Length	3.5 (±0.5)		
		C	Cap Top Width	29 (±1)		
		D	Cap Top Length	18.5 (±1)		
		E	Back of Cap Length	21 (±1)		
UN91251	Green	A	Bottom Width	48 (±2)	100	600
		B	Tie Length	3.5 (±0.5)		
		C	Cap Top Width	29 (±1)		
		D	Cap Top Length	18.5 (±1)		
		E	Back of Cap Length	21 (±1)		
UN91252	Pink	A	Bottom Width	48 (±2)	100	600
		B	Tie Length	3.5 (±0.5)		
		C	Cap Top Width	29 (±1)		
		D	Cap Top Length	18.5 (±1)		
		E	Back of cap length	21 (±1)		
UN91253	White	A	Bottom Width	48 (±2)	100	600
		B	Tie Length	3.5 (±0.5)		
		C	Cap Top Width	29 (±1)		
		D	Cap Top Length	18.5 (±1)		
		E	Back of Cap Length	21 (±1)		

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat, direct sunlight.		
Shelf life	3 Years from the date of manufacturer, if stored as recommended		
Contraindications/Warnings	None		
Markings	Universal Logo, Ref code, Quantity, Product name, Colour, CE Mark, Single use, Latex free, medical device, Non-sterile, keep dry, keep away from sunlight, Barcodes, Recyclable, Lot, Manufacturing & Expiry date, EC-REP, Manufacturer name and address		
Additional information	The product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	26/06/2023	Update to Mediq UK Brand	Ashley Miller
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024