


Product Specifications Sheet

Specification Sheet No	SSN-SCAP(NS)-007		
Product Name	Surgical Caps		
Product Description	Theatre Cap – Elasticated Black (Patterned Material)		
Range	Universal		
Manufacturer	Mediq Healthcare UK Ltd		
Material	Fibre: 69% Viscose, Binder: 31% Polymer based on acrylic ester, acrylonitrile		
Sterilisation	Non Sterile		
Issue Date	21st May 2024	Version No	R02



Intended use	A single use device, intended to be used as a protective barrier when worn on the hands of healthcare providers during examinations/treatment
---------------------	---

Packaging						
Product Code	Colour	Dimensions (cm)			Qty- Inner	Qty – Shipper Carton
UN91230	Blue	A	Full cycle length with stretched elastic	63 (±2)	100	600
		B	Cap height	14 (±0.5)		
		C	Length of Cap Top	20 (±1)		
		D	Width of Cap Top	20 (±1)		
		E	Elastic Length	7-8 (±1)		
UN91231	Green	A	Full Cycle Length with Stretched Elastic	63 (±2)	100	600
		B	Cap Height	14 (±0.5)		
		C	Length of Cap Top	20 (±1)		
		D	Width of Cap Top	20 (±1)		
		E	Elastic Length	7-8 (±1)		
UN91232	Pink	A	Full Cycle Length with Stretched Elastic	63 (±2)	100	600
		B	Cap Height	14 (±0.5)		
		C	Length of Cap Top	20 (±1)		
		D	Width of Cap Top	20 (±1)		
		E	Elastic Length	7-8 (±1)		
UN91233	White	A	Full Cycle Length with Stretched Elastic	63 (±2)	100	600
		B	Cap Height	14 (±0.5)		
		C	Length of Cap Top	20 (±1)		

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Packaging						
Product Code	Colour	Dimensions (cm)			Qty- Inner	Qty – Shipper Carton
UN91234	Yellow	D	Width of Cap Top	20 (±1)	100	600
		E	Elastic Length	7-8 (±1)		
		A	Full Cycle Length with Stretched Elastic	63 (±2)		
		B	Cap Height	14 (±0.5)		
		C	Length of Cap Top	20 (±1)		
		D	Width of Cap Top	20 (±1)		
		E	Elastic Length	7-8 (±1)		

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place. Keep away from heat, direct sunlight.		
Shelf life	3 Years from the date of manufacturer, if stored as recommended		
Contraindications/Warnings	None		
Markings	Universal Logo, Ref code, Quantity, Product name, Colour, CE Mark, Single use, Latex free, medical device, Non-sterile, keep dry, keep away from sunlight, Barcodes, Recyclable, Lot, Manufacturing & Expiry date, Manufacturer name and address		
Additional information	This product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	26/06/2023	Update to Mediq UK Brand replacing DSN 292	Ashley Miller
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024