


Product Specifications Sheet

Specification Sheet No	SSN-SCRUB(NS)-002			
Product Name	Scrub Suits Trousers			
Product Description	Scrub Suits Trousers			
Range	365 Healthcare			
Manufacturer	Mediq Healthcare UK Ltd			
Material	SMS, Blue, with waist tie			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	A non-sterile, two-piece designed to be worn by healthcare staff to facilitate hygienic conditions typically in areas where surgery is performed and is intended to provide a barrier to help prevent cross-contamination of infectious agents between staff and patients. This is a single-use garment.
--------------	--

Packaging						
Product code	Size	Colour	Waist (cm) A	Length (cm) D	Qty- Inner	Qty – Shipper Carton
20-740365/NS	Small	Blue	102±4	104±2	N/A	100
20-741365/NS	Medium		112±4	106±2		
20-742365/NS	Large		116±4	108±2		
20-743365/NS	X Large		120±4	110±2		
20-744365/NS	XX Large		128±4	117±2		

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place. Keep away from heat and direct sunlight.		
Shelf life	5 years from the date of manufacture with the above storage condition		
Contraindications/Warnings	None		
Markings	365 logo, Product name, Quantity, Ref, CE Mark, Latex free, Single use, medical device, Non-sterile, keep dry, keep away from sunlight, Manufacturer name & Address, Barcode, EC-REP, Lot, Manufacturing date & Use by date		
Additional information	The product is latex free and non-sterile.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	26/06/2023	Update to Mediq UK Brand replacing DSN 734	Ashley Miller
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024