


Product Specifications Sheet

Specification Sheet No	SSN-SGO(NS)-003		
Product Name	Surgical Gowns		
Product Description	Standard protection Non-sterile Gown		
Range	365 Healthcare		
Manufacturer	Mediq Healthcare UK Ltd		
Material	Gown: 43g/m ² AR SMS		
Sterilisation	Non-sterile		
Issue Date	21st May 2024	Version No	R02



Remark

- A- Neck circumference
- B- Shoulder width
- C- Cuff length
- D- Sleeve Length
- E- Gown Length
- F- Gown Width

Intended use	Intended to be used to prevent cross contamination between patients and end user
--------------	--

Packaging					
Product code	Size	Dimension (cm)	Colour	Qty- Inner	Qty – Shipper Carton
30-03S365/NS	Large	A = 65 (+2/-0) B = 19 (±1) C = 9.5 (±2) D = 58 (0/-2) E = 115 (±2) F = 156 (±2)	Blue	1 Piece/bag	50 PCS/bag

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	365 Logo, Ref, Size, Product name Image, Quantity, Lot, Barcodes Manufacturing date, Expiry date, Non-Sterile, CE Mark, Single Use, Latex Free, Medical device, EC-REP, Manufacturer Name & Address		
Additional information	This product conforms to standard EN 13795-1:2019 and is latex free.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	30/06/2023	Update to Mediq UK Brand	Racheal begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.