


Product Specifications Sheet

Specification Sheet No	SSN-TA(NS)-001			 A – Applicator Length B – Stick Length C – Stick Width D – Cotton Length E – Cotton Width
Product Name	Applicators			
Product Description	Cotton Tipped Applicator			
Range	Universal			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Cotton tip with Handle Metolose - 60-SH-4000 (as starch) Non-toxic glue			
Sterilisation	Non-Sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	For general use in medical preparation and examinations
--------------	---

Packaging									
Product code	Handle	A (mm)	B (mm)	C (mm)	D (mm)	E (mm)	Items Per Box	Qty- Inner	Qty – Shipper Carton
UN982	White birchwood	150±2	146±1	2.2±0.1	5.5±0.5	14±1	100	20	290
UN984	White Polypropylene	204±2	200±1	3.9±0.1	11-12	35±2	50	10	144
UN986	White or Transparent Polypropylene	76±2	72±0.5	2.5±0.1	5±0.5	12±1	200	24	130

Quality Certification	ISO 13485		
Device Classification	I	Rule	5
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place. Keep away from heat, direct sunlight.		
Shelf life	5 years from the date of manufacture with the above storage condition		
Contraindications/Warnings	None		
Markings	Universal logo, Product name, Quantity, Ref, CE Mark, Single use, medical device, Non-Sterile, keep dry, keep away from sunlight, Do not use if package is damaged, Manufacturer name & Address, Barcode, EC-REP, Lot, Use by date		
Additional information	These products are latex free and non-sterile. Not recommended for use in the ear canal.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	03/07/2023	Update to Mediq UK Brand	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.