


Product Specifications Sheet

Specification Sheet No	SSN-TA(NS)-002			
Product Name	Foam Tipped Applicator/Swab			
Product Description	Foam Swab sticks			
Range	Universal			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Polyurethan head with Polystyrene Handle			
Sterilisation	Non-sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	For general use in medical preparations and examinations
--------------	--

Packaging					
Product code	Length A (mm)	Length B (mm)	Items per Poly Bag	Items Per Pack	Qty – Shipper Carton
UN286	102±2	18±2	N/A	250	4000
UN286-5	102±2	18±2	5pcs	N/A	300

Quality Certification	ISO 13485		
Device Classification	I	Rule	5
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place. Keep away from heat, direct sunlight.		
Shelf life	5 years from the date of manufacture with the recommended storage conditions.		
Contraindications/Warnings	Not recommended for use in the ear canal.		
Markings	Universal logo, Product name, Quantity, Ref, CE Mark, Single use, medical device, Non-Sterile, keep dry, keep away from sunlight, Do not use if package is damaged, Manufacturer name & Address, Barcode, EC-REP, Lot, Use by date		
Additional information	These products are latex free and non-sterile.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	03/07/2023	Update to Mediq UK Brand	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.