


Product Specifications Sheet

Specification Sheet No	SSN-WCD(NS)-001			
Product Name	Wound Care Dressings			
Product Description	Extra Absorbent Dressing Pad			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Non-woven 100% PP, 13±2g/m² Fibre paper, 15±2g/m² Wood pulp, cotton			
Sterilisation	Non-sterile			
Issue Date	21st May 2024	Version No	R02	

Intended use	Used to absorb exudates from wounds. Can be held in place with a retention bandage or adhesive tape.
--------------	--

Packaging					
Product code	Colour	Width ((±1)cm)	Length (cm)	Qty- Inner	Qty – Shipper Carton
1158	White	10	12	100	2,000
1160		20	10	50	1,000
1162		20	20	25	500
1164		20	40	12	240

Quality Certification	ISO 13485		
Device Classification	I	Rule	4
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place away from heat and direct sunlight.		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	None		
Markings	Premier Logo, Ref Code, Quantity, Product name, Dimension, Colour, CE Mark, Single use, Latex free, medical device, Non-sterile, keep dry, keep away from sunlight, Barcodes, Lot, Manufacturing & Expiry date, EC-REP, Manufacturer name and address		
Additional information	The product is latex free		

Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	10/07/2023	Update to Mediq UK Brand replacing DSN 612	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024