

Surgical Drapes

Value & quality, without clinical compromise



Surgical Drapes



Protection & Reassurance

365 Healthcare surgical drapes are trusted by healthcare professionals who require consistency, quality, and confidence in every procedure.

Our drapes are constructed from high-barrier, fluid-resistant materials that help maintain the integrity of the sterile field throughout surgical procedures.

We recognise the critical role that barrier protection plays in preserving a sterile environment, and our products are designed to perform reliably and without compromise.



Fully compliant with EN 13795 standards, our range supports procedures from general surgery through to more complex interventions. Each drape is optimised to enhance procedural efficiency—helping to minimise risk and improve control within the sterile environment.

Our surgical drape range is developed to support infection control, aseptic technique, and overall procedural effectiveness.

Our surgical drapes feature either a 2 or 3-ply laminate construction

These layers are bonded together to provide enhanced strength, absorbency, and repellency, while also reducing the risk of bacterial penetration. Premium adhesives are used to ensure the drapes remain securely in place throughout the procedure.

All 365 Healthcare drapes and sets are blue, a colour chosen to harmonise with the operating theatre environment. This choice minimises glare, improves visibility and contributes to a safe, efficient surgical setting.



Draping Options

The 365 portfolio of surgical drape sets is designed to meet the needs of a wide range of surgical procedures. This versatile range includes the most essential drapes and sets required for many surgical procedures performed regularly.

In addition to these comprehensive sets, many of the drapes are available separately as supplementary options. Popular choices include universal sets, sterile trolley covers, and utility drapes, providing flexibility and convenience for procedural needs.

Benefits

- **High-level strength and absorbency:**

Our surgical drapes are inherently strong due to their material composition they are made of multi-layered, fluid-resistant fabrics and help to safeguard the sterile field.

- **Effective surgical bacterial barrier:**

Provides a robust barrier against bacteria, blocking the passage of microorganisms to keep the surgical site sterile.

- **Lightweight, breathable fabric for comfort:**

Lightweight and breathable fabric ensures comfort for the patient and ease of handling for the surgical team.


- **Procedure-specific draping solutions:**

Tailored to specific surgical procedures, providing optimal coverage and protection.

- **Quality, cost-effective draping solution:**

Offers a reliable and cost-effective solution, balancing quality with budget constraints.

Key Features



Orientation markers
Enhances safety and comfort



Finger Lifts on backing of adhesive
Enhances safety and comfort



Absorbent Panel
When higher fluid levels expected



Set packed in order of use
Supports draping protocols and ease of setting up





Environmental Considerations:

Mediq has shown commitment to sustainability and responsible sourcing by supplying Surgical Drape outer cases which are FSC certified.

The Forest Stewardship Council® (FSC®) is a global, non-profit organisation, providing trusted solutions to help protect the world's forests and tackle today's deforestation, climate, and biodiversity challenges. FSC-certified forests are managed to strict environmental, social and economic standards.


To learn more, visit www.fsc.org

Surgical Draping Options

	Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
General Surgery						
	Universal Set 1	2 x Adhesive edged drapes 2 x Adhesive edged drapes 4 x Hand towels Tape Strip 1 x Zone reinforced trolley cover (sterile field)	150 x 175 cm 75 x 75 cm 30.5 x 48 cm 10 x 50 cm 140 x 190 cm	81-001365	VJD4272	14
	Universal Set 2 with Reinforced Drapes	1 x Reinforced adhesive edged drape (top drape) 1 x Reinforced adhesive edged drape (bottom drape) 2 x Reinforced edged drape (side drape) 4 x Hand towels 1 x Tape strip 1 x Zone reinforced trolley cover (sterile field)	150 x 240 cm 180 x 175 cm 75 x 90 cm 30.5 x 48 cm 10 x 50 cm 140 x 190 cm	81-002365	VJD4273	12
	Universal Set with Reinforced Drape & Mayo	1 x Reinforced adhesive edged bottom drape 1 x Reinforced adhesive edged top drape 2 x Reinforced edged drape 4 x Hand towels 1 x Tape strip 1 x Mayo cover reinforced 1 x Zone reinforced trolley cover (sterile field)	180 x 175 cm 180 x 300 cm 100 x 100 cm 30.5 x 48 cm 10 x 50 cm 80 x 145 cm 140 x 190 cm	81-003365	VJD4274	10
	Incise Drape	1 x Incise Drape	40 x 43 cm Adhesive Area 40 x 35 cm	36588014	ELW414	10
	Incise Drape	1 x Incise Drape	30 x 33 cm Adhesive Area 30 x 25 cm	36588008	ELW875	25




	Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
	Hip Set 1	1 x Reinforced adhesive edged top drape 1 x Reinforced adhesive edged bottom drape 1 x Reinforced split sheet 4 x Hand Towels 1 x Tape strip 1 x Zone reinforced trolley cover (sterile field)	180 x 300 cm 180 x 175 cm 200 x 260 cm 39 x 48 cm 10 x 50 cm 140 x 190 cm	82-001365	VJD4275	8
	Hip Set 3 with Split Sheets & Mayo	1 x Utility drape 1 x Reinforced adhesive edged top drape 1 x Reinforced split sheet 1 x Mayo cover reinforced 4 x Hand towels 2 x Tape strip 1 x Stockinette Large 2 x Zone reinforced trolley cover (sterile field)	75 x 90 cm 180 x 300 cm 200 x 260 (split) 20 x 100 cm 80 x 145 cm 39 x 48 cm 10 x 50 cm 30 x 120 cm 140 x 190 cm	82-003365	VJD4277	8
	Hip Drape with Side Pouches (Elliptical)	1 x Plastic split sheet 1 x Reinforced hip drape (elasticated aperture and tube holders) 4 x Hand towels 2 x Tape strip 1 x Zone reinforced trolley cover (sterile field)	150 x 200 cm (split 20 x 80 cm) 300/230 x 354 cm Aperture 16 x 21 cm 39 x 48 cm 10 x 50 cm 140 x 190 cm	82-004365		8
	Extremity Set 1 with Split sheet	1 x Reinforced extremity drape (elasticated aperture and tube holders) 2 x Hand towels 2 x Tape strip 1 x Plastic split sheet 1 x Zone reinforced trolley cover (sterile field)	200 x 300 cm Aperture 7cm 30.5 x 48 cm 10 x 50 cm 150 x 200 (split 20 x 80 cm) 140 x 190 cm	82-005365	VJD4279	10
	Knee Arthroscopy Drape & Trolley Cover	1 x Arthroscopy drape with fluid collection pouch 1 x Zone reinforced trolley cover (sterile field)	290/230 x 310 cm Aperture 7cm 140 x 190 cm	82-006365	VJD4280	5 sets per inner, 10 sets per outer




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Orthopaedic




DHS Vertical Isolation	Isolation Drape 1 x Fluid collection pouch - large 2 x Fluid collection pouch - Small Incise Film Area Sterile Field Wrap	240 x 330 cm 54 x 49 cm 25 x 29 cm 40 x 30 cm 100 x 100 cm	82-007365	VJD4281	28
------------------------	---	--	-----------	---------	----



T Drape Hand & Foot	1 x Hand drape 1 x CSR wrap (sterile field)	360/160 x 290 cm 100 x 100 cm Aperture 3.5cm	82-008365	VJD4282	16
---------------------	--	--	-----------	---------	----




Foot & Extremity	1 x Foot drape (3.5 cm diameter fenestration) 1 x PSB wrap (sterile field)	200 x 304 cm Aperture 3.5cm 100 x 100 cm	82-011365	VJD4283	18
------------------	---	--	-----------	---------	----



Bilateral Extremity Foot & Leg	1 x Bi-lateral limb sheet (5cm diameter fenestration) 1 x PSB wrap (sterile field)	200 x 305 cm Aperture 5cm 100 x 100 cm	82-014365	VJD4284	18
--------------------------------	---	--	-----------	---------	----




Lower Limb Tourniquet	Tourniquet drape	80 x 45 cm	82-015365/NS	VJD4285	48
-----------------------	------------------	------------	--------------	---------	----




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Orthopaedic




Lower Extremity Foot Drape	1 x Extremity Foot Cover	100 x 35 cm	82-016365	VJD4286	90
----------------------------	--------------------------	-------------	-----------	---------	----




Vertical Isolation Drape - Femoral Nailing	1 x Vertical isolation drape Main drape Incise Area Sterile Field Wrap	240 x 300 cm 25 x 86 cm 100 x 100 cm	82-020365	VJD4250	20
--	---	--	-----------	---------	----



Patient Safety Beach Chair Shoulder	1 x Shoulder drape 1 x PSB wrap (sterile field)	406 x 260 cm Aperture 14 x 13 cm 100 x 100 cm	82-023365	VJD4291	12
-------------------------------------	--	--	-----------	---------	----



Orthopaedic Surgical Bar	1 x Bar drape	280 x 180 cm	82-025365	VJZ153	20
--------------------------	---------------	--------------	-----------	--------	----




Extremity	1 x Reinforced extremity drape (elasticated aperture and tube holders) 1 x PSB wrap (sterile field)	200 x 300 cm Aperture 7 cm 90 x 90 cm	82-027365	VJD4321	20
-----------	--	--	-----------	---------	----



Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Cardiovascular



Lower Vascular	1 x Lower vascular drape	200 x 250 cm	83-004365		20
----------------	--------------------------	--------------	-----------	--	----

Bi-laminate Tri-laminate

PE

PE

Absorbent reinforced panel

Absorbent layer panel

Latex free elasticated material


Adhesive edge

Mesh

Incise area

Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Gynaecology & Obstetrics



Caesarean Section Set


1 x C-section drape with full incise & fluid collection pouch
Incise area
1 x Baby blanket
1 x Baby hat
2 x Velcro strips/tube holders
4 x Hand towels
1 x Zone reinforced trolley cover (sterile field)

257/196 x 325cm
40 x 40cm
110 x 90 cm
12 x 15 cm
2.5 x 30 cm
35.5 x 67 cm
140 x 205 cm

17-006365

VJD4269

5



Caesarean Section Drape (Teardrop Fenestration)


1 x C-Section drape

257/196 x 325 cm
Aperture 22 x 18cm

84-002365

VJD4253

5



Lithotomy Drape Set


2 x Lithotomy leggings
1 x Under buttock drape
1 x Adhesive drape
1 x CSR wrap (sterile field)

75 x 140 cm
75 x 100 cm
90 x 150 cm
75 x 75 cm

84-004365

VJD4254

24



C-Section Drape Set (Teardrop Aperture)

1 x C-Section Drape
2 x Tape Strips
4 x Hand Towels

1 x Table Cover
1 x Baby Blanket
1 x Baby Hat

257/196 x 325 cm
2 x 30 cm
35.5 x 67.5 cm
Aperture
22 x 18cm
140 x 205 cm
90 x 110 cm
12 x 15 cm


84-006365

5




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Urology




TUR Drape set with Trolley Cover	1 x TUR Drape with finger cot and fluid collection pouch 1 x Hand towel 1 x Velcro strip/tube 1 x Zone reinforced trolley cover (sterile field) 1 x Tape Strip	255/155 x 162 cm Aperture 6 cm 30.5 x 48 cm 2.5 x 30 cm 140 x 190 cm 10 x 50 cm	85-001365	VJD4256	14
----------------------------------	--	--	-----------	---------	----

ENT/ Head & Neck




Head & Neck Turban Drape Set	1 x Drape (40 cm long split) 1 x Head drape (75cm long split) 1 x Tape strip 1 x CSR wrap (sterile field)	100 x 105 cm 170/120 x 150 cm 10 x 50 cm 75 x 75 cm	86-002365	VJD4258	40
------------------------------	--	--	-----------	---------	----




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries - Utility Drapes




Utility Drape	1 x Utility drape (2ply)	50 x 60 cm	89-011365	VJD4266	128
---------------	--------------------------	------------	-----------	---------	-----




Utility Drape	1 x Utility drape (2ply)	75 x 90 cm	89-012365	VJD4267	70
---------------	--------------------------	------------	-----------	---------	----



Utility Drape	1 x Utility drape (3ply)	75 x 45 cm	89-014365	VJD4268	94
---------------	--------------------------	------------	-----------	---------	----







Utility Drape	1 x Utility drape (3ply)	75 x 90 cm	89-015365	VJD4293	56
---------------	--------------------------	------------	-----------	---------	----




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Supplementaries - Leggings


	Leggings	Lithotomy Leggings (2 pieces)	75 x 140 cm - 15 cm (folded)	89-001365	ELW926	30
						
	Leggings	Transparent Lithotomy Leggings (2 Pieces)	75 x 120 cm	89-062365	VJD4292	60
						

Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries - Trolley Cover




Trolley Cover	Trolley Cover	100 x 140 cm	89-004365	VJD4260	60
---------------	---------------	--------------	-----------	---------	----




Trolley Cover	Trolley Cover	140 x 140 cm	89-005365	VJD4261	40
---------------	---------------	--------------	-----------	---------	----



Trolley Cover	Trolley Cover	190 x 140 cm	89-006365	VJD4262	38
---------------	---------------	--------------	-----------	---------	----




Trolley Cover	Trolley Cover	210 x 140 cm	89-007365	VJD4263	56
---------------	---------------	--------------	-----------	---------	----




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries - Split Drapes




Split Drape	Plastic Split Sheet	70 x 30 cm (split 20 x 64 cm)	89-031365	VJD4303	56
-------------	---------------------	----------------------------------	-----------	---------	----



Split Drape	Split Sheet with Reinforcement	200 x 260 cm (split 20 x 100 cm)	89-032365	VJD4304	26
-------------	--------------------------------	-------------------------------------	-----------	---------	----




Split Drape	Clear U Drape	120 x 130cm (split 15 x 48 cm)	89-077365	VJZ152	40
-------------	---------------	-----------------------------------	-----------	--------	----




Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries - Adhesive Drapes




Adhesive Drapes	Adhesive Drape	75 x 75 cm	89-019365	VJD4295	72
-----------------	----------------	------------	-----------	---------	----




Adhesive Drapes	Adhesive Drape	75 x 90 cm	89-021365	VJD4296	50
-----------------	----------------	------------	-----------	---------	----



Adhesive Drapes	Adhesive Drape	100 x 100 cm	89-022365	VJD4297	100
-----------------	----------------	--------------	-----------	---------	-----






Adhesive Drapes	Adhesive Drape	150 x 240 cm	89-034365	VJD4305	28
-----------------	----------------	--------------	-----------	---------	----



Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Supplementaries - Adhesive Drapes

	Adhesive Drapes	Adhesive Drape	150 x 240 cm	89-036365	VJD4306	36
	Adhesive Drapes	Adhesive Drape	150 x 175 cm	89-037365	VJD4307	60
	Adhesive Drapes	Adhesive Drape	90 x 90 cm	89-070365	VJD4374	70


Bi-laminate Tri-laminate


PE


PE


Absorbent reinforced panel


Absorbent layer panel


Latex free elasticated material



Adhesive edge


Mesh



Incise area

Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries - Aperture Drapes




Aperture Drapes	Aperture Drape with Adhesive	75 x 90 cm Aperture 6 x 8 cm	89-024365	VJD4298	50
-----------------	------------------------------	---------------------------------	-----------	---------	----



Aperture Drapes	Aperture Drape with Adhesive	50 x 60 cm Aperture 6 x 8 cm	89-025365	VJD4299	180
-----------------	------------------------------	---------------------------------	-----------	---------	-----



Aperture Drapes	Aperture Drape with Adhesive Circular Fenestration	100 x 125 cm Aperture 11cm	89-027365	VJD4300	40
-----------------	--	-------------------------------	-----------	---------	----





Aperture Drapes	Aperture Drape with Adhesive Circular Fenestration	150 x 175cm Aperture 15cm circular	89-051365	VJD4312	50
-----------------	--	---------------------------------------	-----------	---------	----



Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------

Supplementaries - Pouches

	Suction and Diathermy Pouch	Clear Suction and Diathermy Pouch	40 x 35 cm	89-029365	VJD4302	200
	Under Buttock Drape with Fluid Collection Pouch	Under Buttock Drape with Pouch and Fluid Loss Marker	100 x 110cm	89-028365		36

Bi-laminate Tri-laminate

PE

PE

Absorbent reinforced panel

Absorbent layer panel

Latex free elasticated material


Adhesive edge

Mesh

Incise area

Product	Description	Size	Manufacturer Code	NHSSC NPC Code	NHSSC UoM - Box
---------	-------------	------	-------------------	----------------	-----------------


Supplementaries




Drape Accessories	Drape Accessories Tape	10 x 50 cm	89-041365	VJD4310	400
-------------------	------------------------	------------	-----------	---------	-----



Drape Accessories	Stretchable Stockinette Large	32 x 120 cm	89-010365	VJD4265	48
-------------------	-------------------------------	-------------	-----------	---------	----



Mayo Stand Cover	Drape Mayo Reinforced Cover	Cover Dimension 145 x 80 cm Reinforcement Dimension 60 x 75 cm	89-008365	VJD4264	72
------------------	-----------------------------	---	-----------	---------	----



Mediq:

Delivering Excellence in Healthcare Services and Solutions



Product Portfolio

We have a strong and growing brand portfolio across multiple categories.

For example

(this list is not exhaustive)

- Wound care
- Surgical gowns
- Surgical headwear
- Surgical drapes
- Facial protection
- Oxygen & suction tubing
- Procedure packs
- X-Ray detectable swabs



Dedicated account manager & clinical specialists

Account Manager to provide you with support, guidance and advice.

Clinical Specialists that can support care pathways and provide support and advice.

Provide clinical verification of our product portfolio and alternatives.



Implementation excellence & savings analysis

Offering you different options from simple supply to full-service implementation.

Full-service encompasses audit of current products and selection.

Assisting in savings reviews to demonstrate the most cost-effective solutions.



Bespoke Service Level Agreements

We offer adaptable service models, working with you to shape agreements and solutions that are bespoke to your organisation.



at the heart of healthcare

For further information, please contact 01530 830 830 Visit our website www.mediq.co.uk
Browse and order direct from Mediq Online: medbis.mediq.co.uk

©2025 Mediq Healthcare UK Limited. All rights reserved. All other trademarks or brand names are the property of their respective owners. Mediq UK Healthcare Ltd, Mediq House, Arundel Ave, Castle Donington, Derby DE74 2HJ.

Company Registration Number: 00062537, Reference Number: MED36511/0425