



ELIMINATING HAZARDOUS CHEMICAL ACCELERATORS FROM PI SURGICAL GLOVES

Ansell have now eliminated **all chemical accelerators known to cause Type IV allergies or sensitivities** from their PI surgical gloves

Chemical Accelerators	Old Ansell PI	PI-KARE™ Technology
ZDEC (Zinc Diethyldithiocarbamate)	YES	ELIMINATED
ZDBC (Zinc Dibutyldithiocarbamate)	YES	ELIMINATED
SCBC (Sodium Dibutyldithiocarbamate)	YES	ELIMINATED
ZMBT (Zinc Mercaptobenzothiazole)	YES	ELIMINATED
DPTU (Diphenylthiourea)	YES	ELIMINATED
DPG (1,3-Diphenylguanidine)	NO	NO

These substances known to **cause Type IV allergies** or sensitivities are **NOT FOUND** in PI products with **PI-KARE™ Technology**



- Thiurams
- Mercaptobenzothiazole-mix related compounds
- Thiourea related compounds
- Standard or conventional **dithiocarbamates**: ZDEC, ZDBC
- **1,3-diphenylguanidine** (DPG)
- **Cetylpyridinium chloride** (CPC)
- Casein and other animal derivatives
- Phthalates
- (p)-phenylenediamines

Ansell, ® and ™ are trademarks of Ansell Limited or one of its subsidiaries. © 2019 Ansell Limited. All rights reserved.

Know the DIFFERENCE

Learn more about your non-latex options on www.ansellnonlatex.com