


Product Specifications Sheet

Specification Sheet No	SSN-EGL(NS)-004			
Product Name	Examination Gloves (Medical)			
Product Description	White Latex Gloves Powder Free			
Range	Premier			
Manufacturer	Mediq Healthcare UK Ltd			
Material	Natural High Grade Rubber Latex			
Sterilisation	Non-sterile			
Issue Date	1 st Apr 2025	Version No	R03	

Intended use	A single use device, intended to be used as a protective barrier when worn on the hands of healthcare providers during patient examinations/treatment
--------------	---

Packaging									
Product code	Size	Palm Width (mm)	Colour	Palm Thickness (mm)	Length (mm)	Force at Break	Items Per Pack	Qty- Inner	Qty – Shipper Carton
P7900	Extra Small	76 +/- 3	Natural	0.11 +/- 0.02	240	6 N (min)	10	100	1000
P7901	Small	84 +/- 3							
P7902	Medium	94 +/- 3							
P7903	Large	105 +/- 3							

Quality Certification	ISO 13485		
Device Classification	I	Rule	1
Conformity Route	Medical Devices Regulation (EU) 2017/745 Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002)		
Storage conditions	Store in a cool dry place. Keep away from heat and direct sunlight		
Shelf life	5 years from the date of manufacture with the above storage condition		
Contraindications/Warnings	Contains Latex		
Markings	Premier logo, Product name, Quantity, Ref, Nitrile, Colour, Powder free, Size, Lot, Manufacturing date, Use by date, Non-sterile, CE Mark, Single use, contains Latex, Medical device, keep dry, keep away from sunlight, Manufacturer name & Address, Barcode, EC-REP		
Additional information	EN 455 Medical Gloves for single use. Parts 1-4, AQL 1.5 Powder Free, Finger Tip Textured, Ambidextrous Fit.		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	23/06/2023	Update to Mediq UK Brand replacing DSN 249	Racheal Begg
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	21/05/2024	Update the picture and remove old packaging	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024