


Product Specifications Sheet

Specification Sheet No	SSN- FT1358-001		
Product Name	Surgical Gowns		
Product Description	Standard Surgical Gowns		
Range	365 Healthcare		
Manufacturer	Sino Protection (Hefei) Medical Products Co., Ltd		
Material	Gown: Blue AR SMMMMS Wrap: SMS Towel: Polyester scrim		
Sterilisation	EO Sterile		
Issue Date	08 th April 2025	Version No	R03



A = Neck Circumference
B = Sleeve Length
E = Gown Circumference
G = Gown Length

Intended use	Intended to be used to prevent cross contamination between patients and end user
--------------	--

Packaging							
Product code	Size	A (cm) +2/-0	B (cm) +0/-2	E (cm) ±2	G (cm) ±2	Qty- pc/bag	Qty – pc/box
30-02S365	Medium	65	54	145	115	1 x Sterile Gown – Colour Blue, 2 x hand towels, 1 x SMS wrap	30
30-03S365	Large		58	156	122		
30-04S365	Extra Large		62	174	130		
30-05S365	Extra Large, Extra Long		65	176	151		
30-07S365	Extra, Extra Long		58	193	124		

Quality Certification	ISO 13485:2016		
Device Classification	Class Is	Rule	1
Conformity Route	Medical Devices Directive 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 years if stored as recommended		
Contraindications/Warnings	No Known		
Markings	365 Healthcare Logo, Ref, Size, Product name, Quantity, Standard Protection -S, Barcode, NAFDAC number, Lot, Date of Manufacture, Expiry date, Manufacturer name & Address, CE Mark, Single use, Latex free, Medical device, Sterile EO, Do not re-sterilize, Do not use if package is damaged, keep dry, keep away from sunlight, EC-REP, Distributor address		
Additional information	In Compliance with EN 13795-1. The product is latex free		
Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.			

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024



365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that it describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	20/11/2023	New document due to manufacturer's CE and 365 branding	Elizabeth Ogunkolade
R02	21/05/2024	Update Mediq Healthcare UK address	Jack Wong
R03	08/04/2025	Update dimension	Jack Wong

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.

Document Reference	Document Name	Reviewed By	Approved By	Revision Number	Issue Date
DMS 04-004	Product Specification Sheet	Jack Wong	Racheal Begg	R02	03 rd April 2024