



Matrix Park, 1 Western Avenue, Buckshaw Village, Chorley PR7 7NB
T : +44 (0) 1772 299 900 F : +44 (0) 1772 299 901

PRODUCT SPECIFICATION

Product Range		AZO™	
Product Name		AZOWIPE® STERILE	
Product Code	Product Description	Size / Dimensions / Tolerances	Material / Composition / Active Ingredients
81160	AZOWIPE™ STERILE Disinfection wipes for use on hard surfaces (gamma irradiated and triple bagged)	200mm ± 5 mm x 200mm ± 10mm	23gsm 60% Viscose and 40% Polypropylene fibres carded thermo-bonded Active ingredients: 70% v/v isopropanol
Intended Use		Disinfectant wipes for use on hard surfaces	
Classification		Biocide	
Product/Device Single Use		Yes – individual wipes are single use	
Components Latex Free		Yes	
Shelf Life		24 months from date of manufacture	
Storage		Store in a cool, well ventilated area; keep container tightly closed	
Sterilisation Method		Gamma irradiation	
Special Instructions			
Refer to MSDS for safety information			
Pack Qty, Carton Qty, Case Qty, Pack/Pouch content Weight, Case Weight		Pallet Dimensions	
Per pack: 200 wipes Per case: 12 canisters Pack weight: 858g Case weight: 10.3kg		1200 x 1000mm	
Pallet Type	Pallet Configuration	Pallet Weight	
HEAT TREATED STD	Cases per layer: 6 Layers per pallet: 5 Quantity per pallet: 30	310Kg approx	
Packaging Information and Size			
Primary packaging	Secondary Packaging	Tertiary Packaging	
High density white polyethylene canister, PET metal foil seal, white polypropylene lid with a 6mm aperture for wipe dispensing. Canister dimensions 240 x 110mm	Cardboard case Case dimensions: 443 x 335 x 260mm	Grade A HEAT TREATED pallet with corner posts	
Labelling Specification			
Primary (inner):			
Product REF, product name, product description, manufacturer name & address & symbol, expiry date & symbol, LOT number & symbol, storage instructions, pack quantity, active ingredients, barcode, 'do not flush' logo, 'do not macerate' logo, Made in Britain logo, artwork 'IC' reference, Instructions for use.			

It is the responsibility of the user to verify that are using the current released version
HARD COPY UNCONTROLLED UNLESS STAMPED, SIGNED AND DATED BY AUTHORISED PERSONEL



Matrix Park, 1 Western Avenue, Buckshaw Village, Chorley PR7 7NB
T : +44 (0) 1772 299 900 F : +44 (0) 1772 299 901

Secondary (middle):	
Product REF, product name, product description, manufacturer name & address & symbol, expiry date & symbol, LOT number & symbol, keep away from sunlight symbol, keep dry symbol, Barcode containing: (240) product reference (10) LOT number (17) expiry date, "Made in the U.K.", artwork 'IC' reference.	
Tertiary (outer):	
Barcode containing: (240) product reference (10) LOT number (17) expiry date, artwork 'IC' reference, product description, Product REF, Expiry date & symbol. LOT number & symbol, PO number. Total Units Per Box. Total Cases Per Layer, Total Layers Per Pallet, Total Units Per Pallet. Total Cases Per Pallet.	
Standards of Compliance	
BS EN ISO 9001:2008 EU BPR (Regulation 528/2012)	
Approval – Vernacare Limited	
Marketing	
Name DAVID NEVILLE	Signed [Signature]
Position PRODUCT MANAGER	Date 5/9/18
Quality	
Name JONNA HODGKINS	Signed [Signature]
Position LEAD QUALITY ENGINEER	Date 10 Sep 18

It is the responsibility of the user to verify that are using the current released version
HARD COPY UNCONTROLLED UNLESS STAMPED, SIGNED AND DATED BY AUTHORISED PERSONEL