Product Specifications Sheet

Specification Sheet No	SSN-JFRU06-017		
Product Name	Suction		
Product Description	UHS Suction Connecting Tube 3mtr x 7mm FFM		
Range	Universal		
Manufacturer	Romsons International		
Material	Polyvinyl Chloride		
Sterilisation	ETO - sterile		
Issue Date	15th October 2024	**Version No**	R01



Intended use	Urine bag used for collection of urine. Suction Connecting Tube is sterile, single use device, intended for use as an accessory in suctioning, during the surgical procedures

Packaging				
Product Code	Length Approx.	Inside Diameter Approx.	Items Per Pack	Qty – Shipper Carton
UN30027FFM-50	300 cm	7 mm	1	50

Quality Certification	ISO 13485:2016		
Device Classification	Class IIa	**Rule**	4
Conformity Route	Medical Devices Directive 93/42/EEC		
Storage conditions	Store in a cool dry place away from heat and direct sunlight		
Shelf life	5 Years from manufacturing date if stored as recommended		
Contraindications/Warnings	N/A		
Markings	Universal Logo, Ref, Product name, Lot, Barcodes, Expiry date, CE Mark, Do not re-sterilize, Latex free, EO sterile, Medical device, Single use, Do not use if package is damaged, keep dry, keep away from sunlight. DEHP Free, Manufacturer Name & Address, EC-REP, Distributor contact		
Additional information	Product is Latex-free		

Medical Device and Reporting System: Incidents caused by this device must be reported immediately to our offices.

The information provided within this Specification Sheet is correct as of last revision but it is subject to change without notice.

365 Healthcare	Mediq House
T: +44 (0) 1332 413 700	Arundel Avenue, Castle Donington
W: https://mediq.co.uk/	Derby, DE74 2HJ

365 Healthcare is a trading division of Mediq Healthcare UK Ltd
Registered in England No. 62537

Amendment Record

This procedure is reviewed to ensure its continuing relevance to the system and process that is describes. A record of contextual additions or omissions is given below. This revision history lists all changes to this document, the date of each revision, and the reason for each revision.

Revision Number	Revision Date	Reason for Revision	Updated by
R01	15/10/2024	New SSN Mediq Healthcare UK brand Template	Elizabeth Ogunkolade

This document is uncontrolled when printed. The master document is controlled electronically, and document users are responsible for ensuring printed copies are valid prior to use to ensure it represents current practice within Mediq Healthcare UK Ltd.