ONE STEP.
ONE SOLUTION.
FOR 70% IPA DISINFECTION



Why Choose Sani-Cloth® 70?

Sani-Cloth® 70 wipes are fast-drying disinfectant wipes formulated with 70% isopropyl alcohol (IPA) for the rapid disinfection of non-invasive medical devices and non-porous hard surfaces. Sani-Cloth® 70 can also be used to pre-clean or decontaminate medical devices prior to sterilisation.



Test Method	Microorganism	Total Contact Time
EN 16615 EN 13727	Bactericidal	1 min (c) / 2 min (d) 30 sec (c)
EN 16615 EN13624	Yeasticidal	1 min (c) / 2 min (d) 30 sec (c)
EN 14476	Enveloped Viruses	1 min (d)
EN 14348	TB	1 min (c)

(c) Effective disinfection on visibly clean surfaces. For soiled areas, pre-cleaning is recommended to ensure full efficacy.
(d) Proven to clean and disinfect on visibly soiled surfaces – no pre-cleaning needed.

DISINFECTS IN 2 MINUTES

Efficient: Complete disinfection in 2 minutes including Flu in 15 seconds.

PHMB FREE

Free From: Formulated without PHMB, a known skin sensitiser, to minimise the risk of allergic reactions.

SUPPORTS NO GLOVES

Supports Gloves Off: Designed for safe, simple cleaning without the need for gloves.*

MADE IN EUROPE

Minimising our carbon footprint with European manufacturing at our UK & Italian plants.

Sani-Cloth® 70 is compatible with our Point of Use dispensers for improved compliance!



Product Description	Product Code	Wipe Size	Classification	PAO	Shelf life
200ct Canister	E701020	200mm x 200mm	MD Class IIa	3 months	36 months
125ct Canister	E701000	185mm x 133mm	MD Class IIa	3 months	36 months

* For prolonged use and exposure to IPA, wear gloves as per trust policy or if you have a sensitivity to alcohol.

To find out more write to us at: contact@pdi-intl.com

contact@pdi-intl.com +44 (0) 8081 697 945 pdi-intl.com PDI International